# Full Alignment Pyramid Kit

## Instructions for Use

189127

EN	Instructions for Use	2
SL	Navodila za uporabo	12
BG	Инструкции за употреба	22
HR	Upute za uporabu	32
SK	Návod na použitie	42
HU	Használati útmutató	52
EL	Οδηγίες χρήσης	62
LV	Lietošanas pamācība	72
LT	Naudojimo instrukcija	82
ET	Kasutusjuhend	92

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Function.....	6
5 Maintenance .....	6
6 Limitations on Use .....	7
7 Fitting Advice.....	8
8 Technical Data .....	9
9 Ordering Information .....	10

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to the Full Alignment Pyramid Kit.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

## Application

The device is designed for use exclusively as part of a lower limb prosthesis.

Once the device is attached there will be facility for:

- $\pm 6$  mm shift A-P (Anterior - Posterior) / M-L (Medial - Lateral)
- Rotation (up to 360 degrees)
- $\pm 7$  degrees angulation when connected to a Blatchford female pyramid interface.

Constructed in accordance with Blatchford manufacturing recommendations.

Intended for single user only.

## Activity Level

This device is suitable for activity levels 1 to 4; weight limits apply, see *Technical Data*.

Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Contraindications

There are no known contraindications if used in accordance with these instructions.

## Clinical Benefits

- Allows connection of the device to other Blatchford approved prosthetic componentry.

## 2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the device e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Avoid exposure to extreme heat and/or cold.

To reduce the risk of injury due to failure or loosening of the screw connections ensure the bolt threads are cleaned thoroughly before each installation.

The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

Be aware of finger trap hazard at all times.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their practitioner if their condition changes.

Always apply the specified torque values to the screws/ nuts. Never replace the screws with alternative screws. Always use the bolts specified.



Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

# 3 Construction

## Principal Parts

- Dome Aluminum
- Pyramid Titanium
- Bolt Titanium
- Washer Stainless Steel
- Nut Stainless Steel

---

## 4 Function

This device provides correct alignment for the individual user in conjunction with the prosthetic socket. The device allows a secure connection of a Blatchford female pyramid adapter to the flat surface of a Blatchford single bolt system.

---

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider.

Changes in performance may include:

- Instability
- Any unusual noise
- Excessive play or loss alignment

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### **Cleaning**

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

***The remaining instructions in this section are for practitioner use only.***

Maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

Carry out the following routine maintenance at least annually:

- Ensure the nut is tight and secure, if not, remove and clean the treads, apply Loctite and the torque to setting shown.
- Check for visual defects that may affect proper function.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider.

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level, or changes in performance of the device.

---

## 6 Limitations on Use

**This device must only be used with Blatchford modular components.**

### **Intended Life**

A local risk assessment should be carried out based upon activity and usage.

### **Lifting Loads**

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### **Environment**

The product is waterproof to a maximum depth of 1 metre.

This device is suitable for submersion in fresh water only. Light surface corrosion affects neither the function nor the security of this device. However, if it is heavily corroded, stop using and contact your practitioner.

Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts.

Thoroughly rinse with fresh water after use in salt or chlorinated water.

Exclusively for use between -15 °C and 50 °C

(5 °F to 122 °F).

Suitable for submersion

# 7 Fitting Advice

*The instructions in this section are for practitioner use only.*

Avoid extremes of alignment.

Align according to the Instructions for Use supplied with the interfacing componentry.


Adjust the relative positions of components to achieve a comfortable, stable standing position and in conjunction with optimisation of other components' a comfortable, stable gait.

Adjust the alignment (A-P slide, rotation and angulation) in combination with interfacing devices properly accommodate flexion, adduction or abduction. Ensure that the wearer can stand in a stable position with the heel and toe evenly loaded.

### Build Line

The device should be positioned and adjusted to achieve the optimal relative positions (as per manufacturers recommendations) of the various elements of the prosthetic limb build. For optimal/minimal loading the device should be positioned such that the build line is as close as practical to the centre of the device without being at the extremes of angular adjustment.

Symptom	Solution
A recurring noise occurs at the device.	Check the proximal and distal attachments are secure. Tighten according to Construction (section 3) for distal end and in accordance with proximal attachment instructions.
The adapter moves out of position.	User must not use the device until adjusted repaired or replaced.


## 8 Technical Data

Operating and Storage Temperature Range:	-15°C to 50°C (5°F to 122°F)
Component Weight:	103 g
Activity Level:	1–4
Maximum User Weight:	Activity level; 4: 100 kg (220 lb) Activity levels; 1-3: 125 kg (275 lb)
Range of Adjustment:	Angular $\pm 7$ degrees shift dependant on mating component 360° rotational
Attachment Type:	Proximal - Male Pyramid Distal - Blatchford Single Bolt
Build Height:	5 mm

---

### Build Height

## 9 Ordering Information

Part Numbers	
Standard Full Alignment Pyramid Kit	189127

### Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

### CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

### Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

### Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

The above may vary by market; consult your local representative for details.

See the Blatchford website for the current full warranty statement.

## **Reporting of Serious Incidents**

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## **Environmental Aspects**

This product is constructed from recyclable materials. Where possible the components should be recycled in accordance with local waste handling regulations.

## **Retaining the Packaging Label**

The practitioner is advised to keep the packaging label as a record of the device supplied.

## **Trademark Acknowledgements**

Registered trademarks of Blatchford Products Limited.

## **Manufacturer's Registered Address**

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**Blatchford Europe GmbH**
**Endolite India Ltd.**