# Pyramid Laminate Socket Adapter

## Instructions for Use

P189410  
P189410S

EN	Instructions for Use	2
SL	Navodila za uporabo	11
BG	Инструкции за употреба	20
HR	Upute za uporabu	29
SK	Návod na použitie	38
HU	Használati útmutató	47
EL	Οδηγίες χρήσης	56
LV	Lietošanas pamācība	65
LT	Naudojimo instrukcija	74
ET	Kasutusjuhend	83

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Function.....	5
5 Maintenance .....	6
6 Limitations on Use .....	7
7 Fitting Advice.....	8
8 Socket Lamination Procedure.....	8
9 Technical Data .....	9
10 Ordering Information .....	9

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for the practitioner and user unless otherwise stated.

The term device is used throughout to refer to the Pyramid Laminate Socket Adapter.

Please read and ensure you understand these instructions, in particular all safety information and maintenance instructions.

## Application

This device is to be used exclusively as part of a lower limb prosthesis, intended for a single user.

The device is a prosthetic socket adapter which is incorporated into a laminated prosthetic socket during socket manufacture. It offers a Blatchford male pyramid distally for secure attachment to a system that has a Blatchford female pyramid connection.

## Features

- Available in stainless steel (P189410S) or aluminum (P189410)
- A male pyramid adapter and dome, facilitating attachment to a Blatchford pyramid system
- A concave upper surface with peripheral holes and a groove proximal to the dome, facilitating integration into a laminated prosthetic socket

## Activity Level

This device is suitable for Activity Levels 1-4 (weight limits apply, see *Technical Data* section). Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Contraindications

There are no known contraindications if used in accordance with these instructions.

## Clinical Benefit

Provides a connection point from the socket to other prosthetic componentry.

## 2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

P189410 can be used as shower prosthesis but is not designed for prolonged submersion in water. After contact with water wipe dry immediately. If this device comes into contact with salt or chlorinated water, it should be rinsed with fresh water and dried. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

P189410S is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

Be aware of finger trap hazard at all times.

To reduce the risk of injury due to failure or loosening of the bolt connections, ensure the bolt threads are cleaned thoroughly before each installation.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their practitioner if their condition changes.

---

## 3 Construction


### Principal Parts (P189410: three-piece construction)

- Socket Adapter Aluminum
- Male Pyramid Stainless steel
- Screw Steel
- Pyramid Cover Wax

### Principal Parts (P189410S: one-piece construction)

- Socket Adapter Stainless steel
- Pyramid Cover Wax

### Component Identification

---

## 4 Function

The device is designed to be incorporated into the distal end of a custom prosthetic socket using a suitable lamination and resin system. The socket should be shaped by a suitably trained person. Details of the lay up and lamination should take in to account the weight and activity level of the user.

The device allows connection and alignment with the remainder of the system.

---

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, loss of alignment or looseness.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### **Cleaning**

Use a damp cloth and mild soap to clean outside surfaces. DO NOT use aggressive cleansers. Thoroughly dry before use.

***The remaining instructions in this section are for practitioner use only.***

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check the attachment to the Laminated Socket (see *Construction* section).
- Check for defects that could affect proper function.
- Check for corrosion.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider (e.g. significant wear).

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

#### P189410

This device can be used as shower prosthesis but is not designed for prolonged submersion in water. After contact with water wipe dry immediately. If this device comes into contact with salt or chlorinated water, it should be rinsed with fresh water and dried. Light surface corrosion will not affect the function or security of this device. However, if heavy corrosion is evident, stop using this device and contact your practitioner.

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

Suitable for shower use

#### P189410S

This device is waterproof to a maximum depth of 1 meter. Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts. Thoroughly rinse with fresh water after use in salt or chlorinated water.

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

Suitable for submersion

---

## 7 Fitting Advice

*The instructions in this section are for practitioner use only.*

Symptom	Solution
A recurring noise occurs at the device.	Check the attachment between the device and the laminated socket. Confirm all bolts securing the pyramid are securely tightened to the correct torque ratings.
The device moves out of position.	Advise the user they must not use the device until adjusted, repaired or replaced. Inspect device. If components are worn or damaged, then remove and replace with new components. Reposition as necessary and retighten as per <i>Construction</i> section.

---

## 8 Socket Lamination Procedure

*The instructions in this section are for practitioner use only.*


Lay up materials and combinations of materials should be appropriate to the user of the device.

## 9 Technical Data

	<b>P189410</b>	<b>P189410S</b>
Operating and Storage Temperature Range:	-15°C to 50°C (5°F to 122°F)	-15°C to 50°C (5°F to 122°F)
Component Weight:	75 g (2.7 oz)	85 g (3 oz)
Activity Level:	1–4	1–4
Maximum User Weight:	Levels 1–3: 125 kg (275 lb) Level 4: 100 kg (220 lb)	Levels 1–3: 125 kg (275 lb) Level 4: 100 kg (220 lb)
Range of Adjustment:	±7° angular	±7° angular
Proximal Attachment:	Male pyramid (Blatchford)	Male pyramid (Blatchford)
Build Height:	17 mm (35 mm)	17 mm (35 mm)

### Build Height

### Storage and Handling

When storing for prolonged periods, ensure the product is free from moisture and stored at room temperature.

## 10 Ordering Information

Device	Part Number
Pyramid Laminare Socket Adapter (Aluminum)	P189410
Pyramid Laminare Socket Adapter (Stainless Steel)	P189410S

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk).

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a practitioner.

## Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

This product is constructed from recyclable materials. Where possible the components should be recycled in accordance with local waste handling regulations.

## Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

## Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

## Manufacturer's Registered Address

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**Blatchford Europe GmbH**
**Endolite India Ltd.**