# Aluminum Tube

## ∅ 30 mm

### Instructions for Use

903220/200—200 mm long

903220/400—400 mm long

903220/500—500 mm long

EN	Instructions for Use	2
DA	Brugsanvisning	11
NO	Bruksanvisning	20
FI	Käyttöohjeet	29
SV	Bruksanvisning	38
RU	Инструкция протезиста	47
ZH	使用说明	56
AR	إرشادات الاستخدام	65

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Function.....	5
5 Maintenance .....	6
6 Limitations on Use .....	6
7 Fitting Advice.....	7
8 Technical Data .....	8
9 Ordering Information .....	9

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for the practitioner and user unless otherwise stated.

The term *device* is used throughout to refer to the 30 mm diameter Aluminum Tube.

Please read and ensure you understand these instructions, in particular all safety information and maintenance instructions.

## Application

This device is to be used exclusively for lower limb prosthetic fittings, intended for a single user.

The device is used as a connecting element between two Blatchford tube clamp adapters for use in transfemoral and transtibial systems. The device allows length adjustment and, when combined with clamp adapters, rotational adjustment in such systems.

The device is available in three standard lengths (200/400/500 mm) and may be cut further for customized limb alignment and length adjustment to suit the amputee.

## Features

- Lightweight, strong, aluminum construction
- Choice of three standard lengths at 30 mm diameter

## Activity Level

This device is suitable for Activity Levels 1-4 (weight limits apply, see *Technical Data*). Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Contraindications

There are no known contraindications if used in accordance with these instructions.

## Clinical Benefit

- Provides a lightweight connection between prosthetic components.

---

## 2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. excessive play or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Avoid exposure to extreme heat and/or cold. Do not place near any heat source. Do not leave in direct sunshine or inside a car in hot weather.

The device is not designed for prolonged submersion but is suitable for shower immersion and outdoor use. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

The user must not adjust or tamper with the setup of the device.

The user is advised to contact their practitioner if their condition changes.


---

## 3 Construction

### Principal Part

- Tube                      30 mm diameter aluminum tube

### Component Identification

---

## 4 Function

This device provides a connecting element between two Blatchford tube clamp adapters and used together with such adapters, provides rotational and length adjustment in a transfemoral or transtibial system.

Three standard tube lengths are available and are intended to be cut to length by a suitably trained person to suit the requirements of the amputee. This enables a comfortable fit and correct alignment to ensure correct function when using the prosthesis.

---

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, excessive play, loss of alignment or significant wear.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces. DO NOT use aggressive cleansers.

*The remaining instructions in this section are for practitioner use only.*

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check for tube tightness. If loose, remove, clean and re-fit.
- Check for defects that could affect proper function.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider (e.g. significant wear).

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

This device is not designed for prolonged submersion. If this device comes into contact with salt or chlorinated water, it should be rinsed with fresh water and dried. Light surface corrosion will not affect the function or security of this device. However, if heavy corrosion is evident, stop using this device and contact your practitioner.

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).


Suitable for shower use

## 7 Fitting Advice

*The instructions in this section are for practitioner use only.*



### Cutting to Length

The device may be cut to the required length (see *Fitting Length* in Section 8). Ensure the device is cut square and the end has sharp edges removed.

941237

30 mm

Symptom	Solution
A recurring noise occurs at any clamping point.	The device must be secured with no movement. Check that the end of the tube has been cut squarely.
The device moves out of position.	User must not use the device until adjusted, repaired or replaced.

## 8 Technical Data

Operating and Storage Temperature Range: -15 °C to 50 °C  
(5 °F to 122 °F)

Component Weight:  
200 mm long 150 g (5 oz)  
400 mm long 285 g (10 oz)  
500 mm long 360 g (13 oz)

Activity Level: 1–4



Maximum User Weight: 125 kg (275 lb)/1–3  
100 kg (220 lb)/1–4

Attachment Type: Proximal/Distal—ø30 mm tube

Range of Adjustment: 360° axial rotation

Length: 200/400/500 mm

### Fitting Length

### Storage and Handling

When storing for prolonged periods, ensure the product is free from moisture and stored at room temperature.

## 9 Ordering Information

Device	Part Number
30 mm diameter aluminum tube - 200 mm long	903220/200
30 mm diameter aluminum tube - 400 mm long	903220/400
30 mm diameter aluminum tube - 500 mm long	903220/500

### Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

### CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

### Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

### Warranty

This device is warranted for 36 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

### Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

### Environmental Aspects

This product is made from a recyclable material. Where possible, the components should be recycled in accordance with local waste handling regulations.

### Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

### **Trademark Acknowledgements**

Blatchford is a registered trademark of Blatchford Products Limited.

### **Manufacturer's Registered Address**

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**Blatchford Europe GmbH**
**Endolite India Ltd.**