Transfemoral Suction Socket Valve Kit for Flexible Sockets

Instructions for Use

099135

EN	Instructions for Use	2
DA	Brugsanvisning	14
NO	Bruksanvisning	26
FI	Käyttöohjeet	38
SV	Bruksanvisning	50
RU	Инструкция протезиста	62
ZH	使用说明	74
AR	إرشادات الاستخدام	86

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	6
5 Maintenance	6
6 Limitations on Use	7
7 Fitting Advice.....	8
8 Assembly Instructions	9
9 Technical Data	11
10 Ordering Information	12

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to the Transfemoral Suction Socket Valve Kit for Flexible Sockets.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

The device is designed for use exclusively as part of a transfemoral lower limb prosthesis.

A prosthetic valve for use in suction sockets made with flexible thermoplastic socket materials.

Constructed in accordance with Blatchford socket manufacturing recommendations.

Intended for single user only.

Features

- Maintains vacuum within the socket.
- Depression of the central button allows manual ingress or expulsion of air, adjusting the level of vacuum within the socket.
- Facilitates donning/doffing by removal of the central threaded valve.

Activity Level

This device is suitable for activity levels 1 to 4; weight limits apply, see *Technical Data*, Section 9. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

- This device is not suitable for use with rigid socket materials.
- Users must have sufficient balance, dexterity and cognitive function to use the device during donning/doffing and understand how to operate the device.
- Users with poor hand or cognitive function may find donning, operation and cleaning the device difficult.
- Poor hygiene
- This device should not be used in sockets that are not total contact.

Clinical Benefits

- Allows manual air expulsion/ ingress management for different activities and requirements.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the device e.g. restricted movement, non-smooth motion, unusual noises or loss of vacuum should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Avoid exposure to extreme heat and/ or cold.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Be aware of finger trap hazard at all times.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their practitioner if their condition changes.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.


The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

3 Construction

Principal Parts

- Patient Valve Key Nylon or Acetal homopolymer
- Valve Body Assembly Acetal homopolymer
- Nut Acetal homopolymer
- Valve Housing Acetal homopolymer
- Former Acetal homopolymer
- Recessing Disc Chloroprene rubber

Component Identification

4 Function

This device provides a means of maintaining a vacuum within a prosthetic socket. The device allows for manual adjustment of that vacuum and expulsion of air via the central button. The device allows donning (pulling in) of prosthetic sockets via the removal of the central valve.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider.

Changes in performance may include:

- Loss of vacuum
- Unusual noise
- Lack or loss of control/ connection
- Deterioration/ changes to the residual limb
- Tissue damage to the residual limb
- Changes in performance of the device

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

Cleaning the Device

The device is integrated into a prosthetic socket that should be cleaned daily to avoid any build up of bacteria. Take care handling the socket and device to avoid picking up dust, grit and other contamination which could cause irritation to the skin.

Dry thoroughly before use.

Cleaning the Residual Limb

1. Inspect residual limb before and after prosthetic limb use or at least daily. Any deterioration in residuum condition should be reported to the practitioner. If the residual limb has tissue damage, advise your practitioner for advice.
2. Clean the skin daily with an unperfumed, pH balanced soap, rinse the skin with clean water to remove all residues.
3. Dry residual limb.
4. Apply lotion to dry skin, as recommended by a practitioner.
5. Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.

The remaining instructions in this section are for practitioner use only.

Maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually.

- Check that the device functions properly, if not fault find as per *Fitting Advice*.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing the device to corrosive elements such as acids, industrial detergents, bleach or chlorine.

The use of creams or lotions with this device should be used with caution as these can foul the valve.

The product is waterproof to a maximum depth of 1 metre. Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to the device. Thoroughly rinse with fresh water after use in salt or chlorinated water.

Exclusively for use between -15 °C and 50 °C
(5 °F to 122 °F).

Suitable for submersion

7 Fitting Advice

The instructions in this section are for practitioner use only.


Symptom	Remedy
The device is unable to maintain vacuum	<ol style="list-style-type: none">1. Check the valve body is securely located within the prosthetic socket.2. Check the valve housing is seated correctly into the recess.3. Dismantle and clean the valve.

Consider the location of the valve on the definitive socket to suit the individual users unique requirements. Most users having a preference of the valve positioned disto-medially and slightly anteriorly in line with the adductor longus.

8 Assembly Instructions

The instructions in this section are for practitioner use only.


1

The diagram shows a cross-section of a plaster cast with a valve. A 'Latex Cone' is applied to the valve area, and a 'White Nylon Stocking' is pulled over the entire cast. A 'Recessing Disc' is placed on the surface of the cast, secured with double-sided tape. Dashed lines indicate the internal structure of the cast.


1. Locally flatten plaster cast where valve is to be fitted, ensuring area is large enough to accommodate Former 940663.
2. Fit a White Nylon Stocking and Latex Cone over the cast.
3. Apply first vacuum to form Stocking and Latex Cone to cast.
4. Locate Recessing Disc 093503 and retain in position with double sided tape.

2


The diagram shows the 'STAGE 2' assembly. A 'Former 940663' is placed on a piece of 'Double Sided Tape' which is cut around the circumference of the cast. The former is positioned over the valve area, with the tape cut out under the center of the former to ensure the 4 off vacuum holes and the central hole are not covered.

1. Place Former on a length of Double Sided Tape (926401 - 50 mm) and cut around the circumference. Cut out tape under centre of Former to ensure the 4 off vacuum holes and the central hole are not covered.
2. Position Former with double sided tape. **Pin prick** the Latex Cone and Disc around the valve area, through the central hole and the 4 off vacuum
3. Put a sheet of very flexible Socket Material (EVA/low density Polythylene Co-polymer) in oven at 150°C for approximately 13 minutes, turning the sheet once after 6 minutes.

3

1. Drape hot plastic cone over cast and apply second vacuum. Allow to cool.
2. While still on cast centralise a 30 mm dia. Hole Cutter in the impression in the plastic formed by the central hole of the Former.
3. Build up distal end of cast/flexible socket with plaster to the knee interface. Ensure the plaster does not enter hole (plug with plasticine).
4. Pull a latex cone over the cast/flexible socket, followed by a nylon stocking to act as a wick for the vacuum.
5. Fabricate the outer container using a standard drape techniques.
6. On completion remove the outer container and remove the Former from the flexible socket.
7. Insert the Valve Housing 093501 through the hole in the flexible socket. (Use Vaseline if required). This is intended to be a tight fit **do not** enlarge hole. Tighten the Nut 093502 using Locking Ring 940063.
8. **Retighten the valve nut after 24 hours.**

4

1. Using Tool No. 940047 unscrew cap and remove and clean all internal parts prior to reassembly.

9 Technical Data

Part No:	099135 (24 mm dia orifice)
Operating and Storage Temperature Range:	-15°C to 50°C (5°F to 122°F)
Component Weight:	27 g
Activity Level:	1-4

10 Ordering Information

Part Number

099135 (24 mm dia orifice)

Spare Parts	
Patient Valve Key	940046
Valve Body Assembly	099334
Nut	093502
Valve Housing	093501
Former	940663
Recessing Disc	093503

Tools

Description	Part Number
Valve Spanner	940047

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is constructed from recyclable materials. Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.