Female Pyramid Adapter (Short)

Instructions for Use

409041

EN	Instructions for Use	2
DA	Bruksanvisning	11
NO	Bruksanvisning	20
FI	Käyttöohjeet	29
SV	Bruksanvisning	38
RU	Инструкция протезиста	47
ZH	使用说明	56
AR	إرشادات الاستخدام	65

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	6
5 Maintenance	6
6 Limitations on Use	7
7 Fitting Advice.....	8
8 Technical Data	9
9 Ordering Information	9

1 Description and Intended Purpose

These Instructions for Use are intended for the practitioner and user unless otherwise stated.

The term *device* is used throughout to refer to the Female Pyramid Adapter (Short).

Please read and ensure you understand these instructions, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively for modular lower limb prostheses, intended for a single user. The device is designed to accept a Blatchford single-bolt connection such as a Blatchford lower limb housing or Blatchford top housing adapter and allows angular and rotational adjustment between the two limb parts. The single-bolt interface also enables horizontal shift, depending on the type of chassis or housing used.

Features

- Lightweight, strong, stainless steel and aluminum construction
- Blatchford female pyramid
- Blatchford single bolt

Activity Level

This device is suitable for Activity Levels 1-4 (weight limits apply, see *Technical Data*). Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

There are no known contraindications if used in accordance with these instructions.

Clinical Benefit

- Allows connection of prosthetic device with a Blatchford male pyramid to another prosthetic device designed to accept a Blatchford single-bolt connection.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb, e.g. any play, unusual noises or loss of alignment, should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not designed for prolonged submersion or shower immersion but is suitable for outdoor use. Ensure any use of the device complies with the conditions given in *Limitations on Use*.

The user must not adjust or tamper with the setup of the device.

The user is advised to contact their practitioner if their condition changes.

To reduce the risk of injury due to failure or loosening of the bolt connections, ensure the bolt threads are cleaned thoroughly before each installation.

Always apply the specified torque value to the screws.


Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

3 Construction

Principal Parts

- Female Pyramid Adapter Aluminum
- Pyramid Set Screws Plated steel (x4)
- Alignment Bolt Titanium
- Washer Stainless Steel
- Alignment Nut Stainless Steel

Component Identification

4 Function

This device connects a male pyramid and a lower limb chassis or top housing allowing angular and rotational adjustment, plus shift if used with single-bolt interface, between the two limb parts to suit individual limb alignment. When tightened securely, using Loctite and correct torque settings, the device offers a secure way to connect the two parts of the limb.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, movement, significant wear or corrosion that is affecting proper function.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces. DO NOT use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check the set screws for tightness. If loose, remove and clean, then apply Loctite 243 and tighten to the correct torque setting, see *Construction* section.
- Check the alignment nut for tightness. If loose, remove clean, and realign the limb; then apply Loctite 222 and tighten to the correct torque setting, see *Construction* section.
- Check for defects that could affect proper function.
- Check for corrosion.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider (e.g. significant wear).

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).


Suitable for outdoor use

7 Fitting Advice

The instructions in this section are for practitioner use only.

Position as close as practical to the build line.

1. Pass the threaded Alignment Bolt through the single bolt hole in the upper surface of the knee chassis or top housing.
2. Invert the knee and with the Female Pyramid Adapter and bolt supported, pass the washer over the Alignment Bolt followed by the Alignment Nut, ensuring that the raised central lip on the nut engages with the mating recess of the washer.
3. Apply Loctite 222 to the thread of the Alignment Bolt and tighten the Alignment Nut and Bolt to 35 Nm.
4. Place the limb right way up and insert the Male Pyramid from the upper limb into the Female Pyramid Adapter, positioning as required and applying angular and rotational adjustment in order to align the limb.
5. Apply Loctite 243 to the thread of the Set Screws then tighten to 15 Nm.


Symptom	Solution
A recurring noise occurs at the alignment interface.	Check the Alignment Nut is correctly aligned with the washer and is fully tightened. Dismantle and clean the interface. Apply Loctite 222 and tighten to the correct torque setting.
The adapter moves out of position.	Check the Set Screws are fully tightened. Apply Loctite 243 and tighten to the correct torque setting. User must not use the device until adjusted, repaired or replaced.

8 Technical Data

Operating and Storage Temperature Range:	-15°C to 50°C (5°F to 122°F)
Component Weight:	110 g (3.9 oz)
Activity Level:	1–4
Maximum User Weight:	Levels 1–3: 125 kg (275 lb) Level 4: 100 kg (220 lb)
Attachment Type:	Proximal—female pyramid (Blatchford) Distal—3/8 inch UNF nut and bolt
Range of Adjustment:	360° rotation ±7° tilt from vertical
Build Height:	10 mm

Build Height

Storage and Handling

When storing for prolonged periods, ensure the product is free from moisture and stored at room temperature.

9 Ordering Information

Device	Part Number
Female Pyramid Adapter (Short)	409041

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is made from a recyclable material. Where possible, the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.