Female Mini Pyramid Tube Adapter

Instructions for Use

189526 (25 mm) –189727 (30 mm)

| | | |
|----|------------------------|----|
| EN | Instructions for Use | 2 |
| SL | Navodila za uporabo | 12 |
| BG | Инструкции за употреба | 22 |
| HR | Upute za uporabu | 32 |
| SK | Návod na použitie | 42 |
| HU | Használati útmutató | 52 |
| EL | Οδηγίες χρήσης | 62 |
| LV | Lietošanas pamācība | 72 |
| LT | Naudojimo instrukcija | 82 |
| ET | Kasutusjuhend | 92 |

Contents

| | |
|--|----|
| Contents | 2 |
| 1 Description and Intended Purpose | 3 |
| 2 Safety Information | 4 |
| 3 Construction | 5 |
| 4 Function..... | 6 |
| 5 Maintenance | 6 |
| 6 Limitations on Use | 7 |
| 7 Fitting Advice..... | 8 |
| 8 Technical Data | 9 |
| 9 Ordering Information | 10 |

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to the Female Mini Pyramid Tube Adaptor.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively as part of a lower limb prosthesis.

The device connects a Male Mini Pyramid to a 25/30mm diameter tube and allows for angular and rotational adjustment.

Intended for a single user only.

Activity Level

This device is suitable for activity levels 1 to 4; weight limits apply, see *Technical Data* Section 8. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

- There are no known contraindications if used in accordance with these instructions.

Clinical Benefits

- Provides a connection point from a pylon tube to other prosthetic componentry.
- Allows connection of Blatchford Child Pyramid devices.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the device e.g. unusual movement or noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Avoid exposure to extreme heat and/or cold.

To reduce the risk of injury due to failure or loosening of the screw connections ensure the screw threads are cleaned thoroughly before each installation.

The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

Be aware of finger trap hazard at all times.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their practitioner if their condition changes.

Always apply the specified torque value to the screws. Never replace the screws with an alternative screws, always use the screws specified.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

3 Construction


Principal Parts

Part No. 189526 (25 mm)

- Adapter Titanium
- Mini Pyramid Screws Plated steel
- Clamp Screw Plated steel
- Washer Plated steel

Part No. 189727 (30 mm)

- Adapter Titanium
- Mini Pyramid Screws Plated steel
- Clamp Screw Plated steel
- Washer Stainless steel

4 Function

The angular movement and rotation prior to tightening the screws allows for alignment of the limb to suit the individual user.

When tightening securely using Loctite on the pyramid screws only and torque settings shown, the device offers a secure way to connect two parts of the limb.

5 Maintenance

Maintenance must be carried out by competent personnel.

The following routine maintenance is to be carried out at least annually:

- Ensure all screws are secure. If not remove and clean screws, reapply Loctite to the pyramid screws only, tighten all screws to the correct torque settings (see Section 3 *Construction*).
- Check for visual defects that may effect proper function.

The user should be advised that any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Any unusual noises
- Excessive play or loss of alignment

The remaining instructions in this section are for practitioner use only.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check is recommended and signs of wear that may affect function should be reported to their service provider.

Advise the user to inform the practitioner/ service provider of any changes in body weight and/ or activity level and any changes in performance of the device.

6 Limitations on Use

It is recommended that only Blatchford products be used in conjunction with the Female Mini Pyramid Tube adapter, see *Technical Data*. Blatchford Mini Pyramid components must only be used in conjunction with other Blatchford Mini Pyramid components.

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

The product is waterproof to a maximum depth of 1 metre.

Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts.

Thoroughly rinse with fresh water after use in salt or chlorinated water.


Suitable for submersion

Exclusively for use between -15 °C and 50 °C
(5 °F to 122 °F).

7 Fitting Advice

The instructions in this section are for practitioner use only.

Position as close as practical to the build line.


| Symptoms | Remedy |
|--|---|
| <ul style="list-style-type: none"> A recurring noise occurs at the device interfaces. | <ol style="list-style-type: none"> Check security of all alignments screws and tube clamp. Apply loctite as required to set screws (not to the clamping screw) and tighten to the correct torque settings. Ensure mating surfaces are clean and end of tube is cut square. |
| <ul style="list-style-type: none"> The adapter has moved out of position. | <p>Advise the user they must not use the device until adjusted, repaired or replaced.</p> <p>Inspect device. If components are worn or damaged, then remove and replace with new components. Reposition as necessary and retighten to torque settings as shown above, adding loctite where applicable .</p> |

8 Technical Data

| | |
|--|--|
| Operating and Storage Temperature Range: | -15 °C to 50 °C (5 °F to 122 °F) |
| Component Weight: | 75 g (3 oz) |
| Maximum User Weight: | 60 kg (132 lb) |
| Activity Level: | 1-4 |
| Attachment: | Mini Pyramid (Blatchford) |
| Build Height: | 189526 (25 mm) - 12 mm 189727 (30 mm) - 12 mm |
| Range of adjustment: | 360° rotation ±5° angular |
| Connection to: | 189526 - 25 mm Tube 189727 - 30 mm Tube |

Key Dimensions:

9 Ordering Information

| Part Numbers | |
|--|--------|
| Female Mini Pyramid Tube Adaptor 25 mm | 189526 |
| Female Mini Pyramid Tube Adaptor 30 mm | 189727 |

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk.

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is constructed from recyclable materials. Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.