PSPC Cylinder

Instructions for Use

239130

EN	Instructions for Use	2
FR	Instructions d'utilisation	15
DE	Gebrauchsanweisung	28
IT	Istruzioni per l'uso	41
ES	Instrucciones de uso	54
NL	Gebruiksaanwijzing	67
PL	Instrukcje użytkowania	80
PT	Instruções de utilização	93
CS	Návod k použití	106

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	6
5 Maintenance	6
6 Limitations on Use	6
7 Fitting Instructions.....	7
7.1 Upper Attachment.....	7
7.2 ESK+ Demountable Shin	8
7.3 SFESK and ESK/MKL with Anterior Extension Stop	9
7.4 SFESK+ Anterior Extension Stop.....	9
7.5 Length Adjustment	10
8 Adjustment Instructions.....	11
9 Functional Checks.....	11
10 Technical Data	12
11 Ordering Information	13

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to the PSPC (Pneumatic Swing Phase Control) Cylinder.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively as part of a lower limb prosthesis. Intended for a single user.

The device is used to control the swing phase in a transfemoral or knee disarticulation lower limb external prosthesis during normal locomotion. It is to be fitted only as part of a Blatchford modular limb system. (See *Ordering Information*.)

Activity Level

As a general guide this device is recommended for users that have the potential to achieve Activity Level 2 or 3. For further advice refer to the instructions for use of the specific Blatchford limb system the device is used in conjunction with.

Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity level 4 users should adhere to the weight limitations set out in the *Technical Data* section.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

This device may not be suitable for Activity Level 1.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Clinical Benefits

The clinical benefits of the device are determined by the specific Blatchford limb system it is used in conjunction with. (See instructions for use supplied with limb system.)

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

After continuous use the device may become hot to the touch.

Do not place near any heat source. Do not leave in direct sunshine or inside a car in hot weather.

The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk.

Be aware of finger trap hazard at all times

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their practitioner if their condition changes.


Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

3 Construction

Principal Parts

- Piston Rod Titanium
- Top Housing Aluminum
- Bottom Block Aluminum
- Barrel Aluminum
- Silencing Ring Nitrile Rubber

Component Identification

4 Function

The device is a pneumatic cylinder which controls the flexion and extension resistance of a limb during swing phase. The device is fitted with two manually adjustable needle valves to allow customization of resistances in both directions.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or restricted/excessive rotation or significant wear.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check smoothness of operation.
- If necessary, relubricate piston rod surface using a drop of STP oil treatment.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider.

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Suitable for outdoor use

Exclusively for use between -10 °C and 50 °C (14 °F to 122 °F).

7 Fitting Instructions

The instructions in this section are for practitioner use only.

Use appropriate health and safety equipment at all times including extraction facilities.



Be aware of finger trap hazard at all times.

Fitting Kit: ESK+ Demountable Shin 019758

7.1 Upper Attachment


Fitting Trunnion Assembly

7.2 ESK+ Demountable Shin

The instructions in this section are for practitioner use only.

 Remove silencing ring and discard prior to fitting cylinder into either ESK+ shin.


 **Caution:** On assembly, ensure the shoulders of the lower pivots align correctly with the shin bearings.

ESK+ Stanceflex shown as example.

7.3 SFESK and ESK/MKL with Anterior Extension Stop

The instructions in this section are for practitioner use only.

When used with SFESK or ESK/MKL with Anterior Extension Stop, the stop should be adjusted as follows:

1. Slacken Trunnion Clamping Screw.
2. Using a 6 mm Spanner rotate piston rod until Knee Lock Patella locks into place without rock and with ease.
3. Retighten the Trunnion Clamping Screw.
4. Remove Patella (except MKL).
5. Using a strip of paper 5 mm wide, place on top of the front buffer. Screw down front extension stop until paper can be pulled out with a slight resistance. A light extension movement should be applied to the top of the socket whilst carrying out this adjustment.
6. Replace the Patella (except MKL).






7.4 SFESK+ Anterior Extension Stop

When fitted with a PSPC, the ESK+ extension stop should be wound fully home. No adjustment is necessary.

7.5 Length Adjustment


The instructions in this section are for practitioner use only.

For modules supplied with control cylinder already assembled into the shin, the following has already been carried out.

8 Adjustment Instructions

The instructions in this section are for practitioner use only.

Setting The Flexion And Extension Valve

The Flexion and Extension Valves should be adjusted to suit the individual wearers gait characteristics. Turn each adjuster within the range specified below, clockwise to increase resistance and anti-clockwise to reduce resistance.

To ensure the limb extends fully always start of with the Extension resistance at a minimal value.

! User setup should initially be done with the support of parallel bars.

Approximate Range of Adjustment

	Minimum	Maximum
Extension	3 Rings Showing	No Rings Showing
Flexion	3 Rings Showing	No Rings Showing

9 Functional Checks

The following manual checks may be made to ensure correct function of the cylinder non return valves and seals.


Valve Positions	Outstroke	Instroke
Extension Resistance Minimum Flexion Resistance Maximum	Free	Cushioned
Extension Resistance Maximum Flexion Resistance Minimum	Cushioned	Free

10 Technical Data

Operating and Storage Temperature Range:	-10 °C and 50 °C (14 °F to 122 °F)
Component Weight:	252 g
Activity Level:	2-3 (4)
Maximum User Weight:	
Activity Levels 2-3	125 kg (275 lb)
Activity Level 4	100 kg (220 lb)
Build Height:	See Diagrams Below.

Fitting Length

Suitable for 160 mm x 30 mm geometry knees.

11 Ordering Information

Compatible Blatchford Knee Systems

ESK+ Knee Module	Part Number
Supplied without PSPC Cylinder Fitted	ESKSF30
	ESKSFPPYR
	ESK4B30
	ESK4BPYR
Supplied with PSPC Cylinder Fitted	ESKSF30
	ESKSFPPYR
	ESK4BP30
	ESK4BPPYR
	ESKSBMKLP30
	ESKSBMKLPPYR

Fitting Kit

PSPC	019758
------	--------

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

The above may vary by market; consult your local representative for details.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

To help prevent potential harm to the environment or to human health from uncontrolled waste disposal, Blatchford offers a take-back service. Please contact Customer Services for details.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.