# VTM

## Instructions for Use

VTM1 – VTM4

EN	Instructions for Use	2
SL	Navodila za uporabo	13
BG	Инструкции за употреба	24
HR	Upute za uporabu	35
SK	Návod na použitie	46
HU	Használati útmutató	57
EL	Οδηγίες χρήσης	68
LV	Lietošanas pamācība	79
LT	Naudojimo instrukcija	90
ET	Kasutusjuhend	101

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Function.....	5
5 Maintenance .....	6
6 Limitations on Use .....	6
7 Bench Alignment.....	7
7.1 Static Alignment.....	7
7.2 Dynamic Alignment .....	7
8 Axial Spring Replacement.....	9
9 Technical Data .....	10
10 Ordering Information .....	11

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to VTM (Vector Thrust Male).

**Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.**

## Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

## Activity Level

This device is recommended for users that have the potential to achieve Activity Level 3 or 4. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

## Contraindications

This device might not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding maintenance.

## Clinical Benefits

- Reduced loading rate on the residual limb, compared to rigid pylons
- Reduced back pain when twisting, compared to rigid pylons
- Increased walking speed, compared to rigid pylons

# Axial Spring Selection Chart

		User Weight									
Impact	lb	100-115	116-130	131-150	151-170	171-195	196-220	221-255	256-285	286-325	326-365
	kg	44-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147	148-166
LOW		1	1	1	2	2	3	3	4	4	4
MODERATE		1	1	2	2	3	3	4	4	4	
HIGH		1	2	2	3	3	4	4	4		
Axial spring rating											

## Impact levels:

- Low** Daily walking and occasional sports such as golf and hiking.
- Moderate** Aggressive walking, frequent or daily sports such as jogging.
- High** Daily activities such as distance running, climbing, lifting and carrying heavy objects for vocational purposes.
- Extreme** Rigorous daily walking at high cadence, competitive athlete in track and field type events, has vocation that requires jumping or loading of the prosthesis beyond normal levels.  
(Extreme users should contact our sales team.)

## 2 Safety Information

**This warning symbol highlights important safety information which must be followed carefully.**

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

The user should contact their clinician if their condition changes.


Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

# 3 Construction

## Principal Parts:

- Housing (Aluminium Alloy)
- Pyramid Tube (Titanium)
- Axial Spring (Titanium)
- Screws (Titanium/Stainless Steel/Nylon)
- Bearings (Reinforced Polymer)

## Component Identification

# 4 Function

The device comprises a housing assembly which contains a spring component that allows the Pyramid Tube to deflect and/or rotate through a limited range under an axial and or torsional load. When the load is removed the spring returns the Pyramid Tube to its original position. The effect of this movement is to attenuate any shock forces that would otherwise be transmitted to the stump socket interface both in axial and rotational directions.

---

## 5 Maintenance

The wearer should be advised:

- Any changes in performance of this device must be reported to the practitioner  
Changes in performance may include:
  - Increased or reduced axial or torsional stiffness
  - Excessive play in bearings
  - Any unusual noises

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces. Do not use aggressive cleansers.

*The remaining instructions in this section are for practitioner use only.*

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

It is recommended that annually the device is checked for visual defects that may affect proper function and that the bearing surface is re-lubricated. Use 'Sapphire Endure' (928015) or equivalent PTFE loaded grease and follow the instructions in "**Bench Alignment**" section for dismantling and re-assembly.

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

Avoid exposing the device to corrosive elements such as water, acids and other liquids.

Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.


Suitable for outdoor use

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

## 7 Bench Alignment

*The instructions in this section are for practitioner use only.*

### 7.1 Static Alignment

#### Set-Up Length

#### Note:

**Shin blade should be vertical when used with a 10 mm heel height shoe.**

Lower heel height shoes may require the use of a packing piece below the heel.

Cut shin to approximate length, allowing 10-15 mm for final trimming. Note trial length will typically include 10 mm to allow for axial compression of heel, toe and spring (axial).

Use a VT Cutting and Drilling Jig (Part No. 941256) to fit the device to an Elite Blade.

#### Build Line

This should fall 1/3 of the foot length from the heel (AP view shown).

### 7.2 Dynamic Alignment

#### Coronal Plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

#### Sagittal Plane

Check for smooth transition from heel strike to toe off. Ensure also that when standing the heel and toe are evenly loaded and that both are touching the floor.

#### Transfemoral Alignment

Align trans-femoral devices according to fitting instructions supplied with the knee, keeping the build line relative to the device as shown.

**Fitting Advice:**

The springs of the device are designed to work together with the heel and toe springs of the Elite Blade to give smooth progression for most users.

The axial spring controls both compressive and torsional resilience. Typical axial movement should be between 3-6 mm for normal walking. Alternative stiffness springs are available if required.

The device is designed to be interchangeable with existing Elite Blade Top Housings that have been fitted with the shifting Male Pyramid option (Part No. 539501 or 539601). No further length adjustment should be necessary.

For additional fitting advice see Elite Blade Instructions for Use (Part No. 938295S).

# 8 Axial Spring Replacement

*The instructions in this section are for practitioner use only.*

Remove the device from the Elite Blade. Lubrication of the bearing surface may be carried out using 'Sapphire Endure' (928015) or equivalent PTFE loaded grease.

**1** Remove the rotation stop

**2** Undo and remove Pyramid screw and withdraw pyramid tube

**5**

Note:  
It may be necessary to partly unscrew the Pyramid Screw as shown and gently tap the screw head to disengage the anti-rotation pins.

**3** Using a long reach hex key undo the retaining screw and remove axial spring

**4**

axial spring

retaining screw

spring rate indicator

**4** Assembly is the reverse of the above, taking care to orientate the Pyramid Tube aperture to ensure access for the rotation stop. Mark the bottom of the carrier with the appropriate rating of spring [ 1-4 ] using a permanent marker (see table).

Axial Spring rating indicator


The  $\beta$  symbol may be used to indicate the chosen foot spring category

Indicator	Spring Rate	Spring (Part No.)
	1	539058
	2	539059
	3	539060
	4	539061

## 9 Technical Data

Operating and Storage Temperature Range:	-10 °C to 50 °C (14 °F to 122 °F)
Component Weight (size 26):	450 g (1 lb)
Activity Level:	3–4
Maximum User Weight:	166 kg (365 lb)
Proximal Alignment Attachment:	Male Pyramid (Blatchford)
Range of Adjustment:	±7° Angular
Maximum Rotation	± 15°
Typical Axial Movement	3-6 mm [max. 10 mm]

### Fitting Length:

# 10 Ordering Information

Product Code	Fitted Spring Rate
VTM1	1
VTM2	2
VTM3	3
VTM4	4

Item	Part Number
Spring - Rate 1	539058
Spring - Rate 2	539059
Spring - Rate 3	539060
Spring - Rate 4	539061
Sapphire Endure	928015

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## Warranty

This device is warranted for 36 months - foot shell 12 months - glide sock 3 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licences and exemptions.

See Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

## Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

## Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH.

## Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

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**Blatchford Europe GmbH**
**Endolite India Ltd.**