SuperSACH

Instructions for Use

P529101–P529114

EN	Instructions for Use	2
SL	Navodila za uporabo	13
BG	Инструкции за употреба	24
HR	Upute za uporabu	35
SK	Návod na použitie	46
HU	Használati útmutató	57
EL	Οδηγίες χρήσης	68
LV	Lietošanas pamācība	79
LT	Naudojimo instrukcija	90
ET	Kasutusjuhend	101

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	6
5 Maintenance	6
6 Limitations on Use	7
7 Bench Alignment.....	8
7.1 Static Alignment.....	8
7.2 Dynamic Alignment	8
7.3 Transfemoral Alignment.....	8
8 Fitting Advice.....	9
8.1 Using the Adjuster Screw	9
8.2 Bonding a Cosmesis	9
9 Technical Data	10
10 Ordering Information	11

1 Description and Intended Purpose

These Instructions for Use are intended for the practitioner and user unless otherwise stated.

The term *device* is used throughout to refer to the SuperSACH.

Please read and ensure you understand these instructions, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively as part of a lower limb prostheses with an interfacing Blatchford product, intended for a single user.

This device provides the user with a cushioned heel and a flexible keel. The adjustable heel screw allows heel cushion adjustment for an individual fit.

Features

- Lightweight polyurethane foam foot shell
- Integrated pyramid insert
- Adjustable heel cushion
- Sandal toe

Activity Level

This device is recommended for Activity Level 1 or 2 users (weight limits apply, see *Technical Data*). Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Contraindications

This device may not be suitable for competitive sports events, as these types of user may be better served by a specially designed prosthesis optimized for their needs.

Clinical Benefits

- Low heel stiffness at weight acceptance gives good stability for lower mobility patients
- Long-lasting for low mobility users
- High satisfaction rates amongst bilateral users

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted or excessive movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, competitive running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling or running is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner that has attended an approved training course.

The user should be advised to contact their practitioner if their condition changes.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

To minimize the risk of slipping and tripping, appropriate footwear that fits securely onto the foot shell must be used at all times.


The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

3 Construction

Principal Parts

- Foot Shell Polyurethane (PU) foam
- Adjuster Screw (M10x40) Nylon
- Pyramid Insert Stainless steel
- Keel Nylon ST
- Dome Aluminum alloy

Component Identification

4 Function

The device comprises a lightweight thermoplastic keel which gives a gentle rolling action through stance. The keel lies within a standard closed foot shell profile and features a sandal toe for a more natural appearance.

An adjuster screw in the heel allows the user to customize the hardness of the heel cushion for individual requirements.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness, changes in heel cushioning, restricted/excessive movement, significant wear or excessive discoloration from long term exposure to UV.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces. DO NOT use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check for damage or wear and replace if necessary.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level

If this device is used excessively, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

This device is waterproof to a maximum depth of 1 meter. Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear. Thoroughly rinse with fresh water after use in salt or chlorinated water.

Foot products must be adequately finished to prevent water ingress into the foot shell where possible. If water enters the foot shell, the limb should be inverted and dried before further use.


Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

Suitable for submersion

7 Bench Alignment

The instructions in this section are for practitioner use only.

7.1 Static Alignment

With the flexion accommodated the build line should pass through the center of the socket in the sagittal (A-P) plane and through the SuperSACH foot, as shown. Make sure you accommodate the heel height of the shoe.

Build Line

This should fall 1/3 of the foot length from the heel.

7.2 Dynamic Alignment

Coronal Plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

Sagittal Plane

Check for smooth transition from heel strike to toe-off. Ensure also that when standing the heel and toe are evenly loaded and that both are touching the floor.

7.3 Transfemoral Alignment

Align transfemoral components according to fitting instructions supplied with the knee, keeping the build line relative to this device, as shown.

* Allow for user's own footwear

8 Fitting Advice

The instructions in this section are for practitioner use only.

Symptoms	Reason	Remedy
<ul style="list-style-type: none">• Sinking at heel strike	Heel too soft	Stiffen the heel cushion
<ul style="list-style-type: none">• Rapid transition from heel strike through stance phase• Difficulty controlling heel action (e.g. foot jars into mid-stance)• Foot feels too rigid	Heel too hard	Soften the heel cushion

8.1 Using the Adjuster Screw


The adjuster screw is used to adjust the stiffness of the heel cushion.

Note... Users who weigh between 100kg and 125kg may need maximum heel stiffness. Users who weigh 60kg or less may need minimum heel stiffness.

To stiffen the heel cushion:

- Turn the adjuster screw clockwise.


Note... The heel cushion is set to maximum stiffness when the screw is fully engaged.

To soften the heel cushion:

- Turn the adjuster screw counterclockwise.

Note... The heel cushion is set to minimum stiffness after 4 counterclockwise turns from the maximum stiffness setting.

8.2 Bonding a Cosmesis

Before you begin, lightly abrade and clean the mating surfaces with a suitable solvent cleaner.


938097

For instructions on fitting a cosmesis

Use a contact adhesive (Evo Stik 528, Thixofix or an equivalent)


Use appropriate health and safety equipment at all times including extraction facilities

9 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight (size 26):	449 g (15.1 oz)
Activity Level:	1–2
Maximum User Weight:	125 kg (275 lb)
Attachment Type:	Blatchford Male Pyramid Adapter
Range of Adjustment:	$\pm 7^\circ$ tilt from vertical
Build Height:	70–80 mm
Heel Height:	19 mm

Build Height

Storage and Handling

When storing for prolonged periods, ensure the product is free from moisture and stored at room temperature out of direct sunlight (protected from UV).

10 Ordering Information

SuperSACH

Size (cm)	Part number — left	Part number — right
22	P529101	P529102
23	P529103	P529104
24	P529105	P529106
25	P529107	P529108
26	P529109	P529110
27	P529111	P529112
28	P529113	P529114

Ankle Cosmesis

Size (cm)	Part number — left	Part number — right
All sizes	P561081	P561081

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is made from a recyclable material. Where possible, the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.