# Full Alignment Kit

## Instructions for Use

189806, 189805, 189807  
189906

EN	Instructions for Use	2
DA	Brugsanvisning	12
NO	Bruksanvisning	22
FI	Käyttöohjeet	32
SV	Bruksanvisning	42
RU	Инструкция протезиста	52
ZH	使用说明	62
AR	إرشادات الاستخدام	72

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Function.....	6
5 Maintenance .....	6
6 Limitations on Use .....	7
7 Fitting Advice.....	8
8 Technical Data .....	9
9 Ordering Information .....	10

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for the practitioner and user unless otherwise stated.

The term *device* is used throughout to refer to the Full Alignment Kit

Please read and ensure you understand these instructions, in particular all safety information and maintenance instructions.

## Application

This device is to be used exclusively for modular lower limb prostheses, intended for a single user.

This device is designed to connect the socket dome to a knee chassis. 189805, 189806 and 189807 are for above knee (AK) only. 189906 is for below knee (BK) only. Sockets must have Blatchford single bolt alignment. The device allows optimal alignment to be achieved when using a socket manufactured with single bolt construction.

## Features

- Lightweight, strong, aluminum/titanium construction

## Activity Level

This device is suitable for Activity Levels 1-4 (weight limits apply, see *Technical Data*). Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Contraindications

There are no known contraindications if used in accordance with these instructions.

## Clinical Benefit

- Aids in setting correct prosthetic alignment during prosthetic socket manufacture.

---

## 2 Safety Information

**This warning symbol highlights important safety information which must be followed carefully.**

Any changes in the performance or function of the limb, e.g. excessive play or unusual noises, should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

The user must not adjust or tamper with the setup of the device.

The user is advised to contact their practitioner if their condition changes.

To reduce the risk of injury due to failure or loosening of the bolt connections, ensure the bolt threads are cleaned thoroughly before each installation.


Always use Loctite and apply the specified torque value to the bolts. Never use an alternative bolt.

# 3 Construction

## Principal Parts

- Socket Plate/Bolt Assembly Aluminum/titanium
- Friction Disk Synthetic rubber
- Upper Alignment Coupling Aluminum
- Spacer Polyarylamide
- Washer Stainless steel
- Alignment Nut Stainless steel

## Component Identification

---

## 4 Function

This device connects a single bolt socket to the knee chassis, providing secure attachment.

109806, 189805 and 189807 are used for transfemoral prosthesis applications.

189906 is used for transtibial prosthesis applications.

---

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises or significant wear. There should be no movement at this interface.

Inform the practitioner/service provider of any changes in body weight and/or activity level. Look for and report any signs of corrosion.

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces. DO NOT use aggressive cleansers.

***The remaining instructions in this section are for practitioner use only.***

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check the alignment nut for tightness. If loose, remove, clean, and realign limb; then apply Loctite 222 and tighten to the correct torque setting, see *Construction* section.
- Check for defects that could affect proper function.
- Check for corrosion.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

This device is waterproof to a maximum depth of 1 meter. Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts. Thoroughly rinse with fresh water after use in salt or chlorinated water.

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).


Suitable for submersion

## 7 Fitting Advice

*The instructions in this section are for practitioner use only.*

Position as close as possible to the build line.

1. With the socket plate inside the socket, pass the threaded Alignment Bolt through the hole in the end of the socket.
2. On the outside of the socket, pass the Alignment Coupling and Friction Disk over the bolt followed by any spacers used, ensuring that the friction surface of the spacer is distal and touches the upper surface of the knee.
3. Pass the Alignment Bolt through the centre of the single bolt hole in the upper surface of the knee.
4. Invert the knee and with the socket plate and bolt supported, pass the washer over the central bolt followed by the nut, ensuring that the raised central lip on the nut engages with the mating surface of the washer.
5. Apply Loctite 222 to the thread of the Alignment Bolt. Adjust and align the limb according to the IFUs supplied with the knee, then tighten the Alignment Nut and Bolt to 35 Nm.
6. Seal with silicone for waterproofing.

Symptom	Solution
A recurring noise occurs between the socket and the knee.	Ensure that the socket interface has been correctly formed and is the correct thickness. The Alignment Nut must be tightened. Apply Loctite and tighten to the correct torque setting.
The adapter moves out of position.	User must not use the device until adjusted, repaired or replaced.


## 8 Technical Data

	<b>189806, 189805, 189807</b>	<b>189906</b>
Operating and Storage Temperature Range:	-15°C to 50°C (5°F to 122°F)	-15°C to 50°C (5°F to 122°F)
Component Weight:	183 g (6.5 oz) 61 g (2.2 oz) 93 g (3.3 oz)	122 g (4.3 oz)
Activity Level:	1–4	1–4
Maximum User Weight:	125 kg (275 lb)/1–3 100 kg (220 lb)/1–4	125 kg (275 lb)/1–3 100 kg (220 lb)/1–4
Attachment Type:	Proximal—socket Distal—3/8 inch UNF nut and bolt	Proximal—socket Distal—3/8 inch UNF nut and bolt
Range of Adjustment:	360° axial rotation ±6° tilt from vertical ±5 mm shift A-P or M-L	360° axial rotation ±6° tilt from vertical ±5 mm shift A-P or M-L
Build Height:	20 mm (standard) 10 mm 30 mm	20 mm

### Fitting Length

The overall height of the device is either 10, 20 or 30 mm, depending on the part to be fitted. Appropriate allowance must be taken into account in the build height of the limb. There are no height adjustment facilities available with this device.

### Storage and Handling

When storing for prolonged periods, ensure the product is free from moisture and stored at room temperature.

---

## 9 Ordering Information

Device	Part Number
Full Alignment Kit TF (20 mm - standard)	189806
Full Alignment Kit TF (10 mm)	189805
Full Alignment Kit TF (30 mm)	189807
Full Alignment Kit TT (20 mm)	189906

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

This product is made from a recyclable material. Where possible, the components should be recycled in accordance with local waste handling regulations.

## Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

## Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

## Manufacturer's Registered Address

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**Blatchford Europe GmbH**
**Endolite India Ltd.**