Multiflex Pyramid Ankle

Instructions for Use

409353-409355

| | | |
|----|----------------------------|----|
| EN | Instructions for Use | 2 |
| FR | Instructions d'utilisation | 12 |
| DE | Gebrauchsanweisung | 22 |
| IT | Istruzioni per l'uso | 32 |
| ES | Instrucciones de uso | 42 |
| NL | Gebruiksaanwijzing | 52 |
| PL | Instrukcja użytkowania | 62 |
| PT | Instruções de utilização | 72 |
| CS | Návod k použití | 82 |

Contents

| | |
|--|----|
| Contents | 2 |
| 1 Description and Intended Purpose | 3 |
| 2 Safety Information | 4 |
| 3 Construction | 5 |
| 4 Function..... | 5 |
| 5 Maintenance | 6 |
| 6 Limitations on Use | 6 |
| 7 Bench Alignment..... | 7 |
| 7.1 Additional Components | 7 |
| 7.2 Static Alignment..... | 7 |
| 7.3 Dynamic Alignment | 7 |
| 8 Fitting Advice..... | 8 |
| 9 Technical Data | 9 |
| 10 Ordering Information | 10 |

1 Description and Intended Purpose

These Instructions for Use are intended for the practitioner and user unless otherwise stated. The term *device* is used throughout these instructions for use to refer to Multiflex Pyramid Ankle. Please read and ensure you understand these instructions, in particular all safety information and maintenance instructions.

Application

This device is an ankle unit to be used exclusively as part of a lower limb prosthesis, intended for a single user.

The device provides the user with comfortable multi-axial motion, which enhances stability when used in combination with a Multiflex foot. Changeable ankle ball and snubbers provide the ability to customize the ankle according to the user's weight and Activity Level.

Activity Level

This device is suitable for Activity Level 1–2 (weight limits apply, see *Technical Data* section).

As an exception to our recommendation, Activity Level 3 and Activity Level 4* users may only use this device based upon the outcome of a prosthetic assessment.

*Maximum user weight: 100 kg (220 lb).

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as, curbs, stairs or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

This device is not suitable for competitive sports events. This type of user will be better served by a specially designed prosthesis optimized for their needs.

Clinical Benefits

- Low stiffness at weight acceptance leads to early foot-flat and greater stability for lower mobility users, compared to single axis feet and Solid-Ankle-Cushioned-Heel feet
- No loss of stability for bilateral users, compared to rigid ankles
- Easier to walk on uneven ground for bilateral users, compared to rigid ankles
- Easier to walk on slopes for bilateral users, compared to rigid ankles
- Easier to walk longer distances for bilateral users, compared to rigid ankles
- Improved stance phase timing symmetry, compared to rigid ankles
- Reduced sound limb loading, compared to rigid ankles

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted or excessive movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Any excessive changes in heel height after set up will adversely affect limb function and should be immediately reported to your service provider.

Do not place near any heat source. Do not leave in direct sunshine or inside a car in hot weather.

The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately. Ensure any use of the device complies with the conditions given in *Limitations on Use*.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Be aware of finger trap hazard at all times

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their clinician if their condition changes.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.


Make sure that no liquids enter the device during donning or normal use.

3 Construction

Principal Parts

- Pyramid Insert Stainless Steel
- Ankle Housing Aluminum
- Ankle Ball Aluminum/Natural Rubber
- Collet Aluminum
- Internal Circlip Plated Carbon Steel

Component Identification

4 Function

The device comprises an ankle housing with a rubber ball joint. When used with an appropriate snubber kit, these components provide the user with a resilient multi-axial ankle motion for secure foot placement at heel strike and a smooth progression through the gait cycle. The ball joint primarily controls plantarflexion and the snubber controls dorsiflexion support.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or increased range of motion.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check the ankle for excessive play
- Check the tightness of the socket cap screw

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider.

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.


Exclusively for use between -10°C and 50°C (14°F to 122°F).

Suitable for outdoor use

7 Bench Alignment

The instructions in this section are for practitioner use only.

7.1 Additional Components

Follow the alignment instructions that are supplied with any additional components.

7.2 Static Alignment

With the flexion accommodated the build line should pass through the center of the socket in the sagittal (A-P) plane and through the center of Multiflex Pyramid Ankle. Heel height of the shoe should be accommodated and adjusted at the serrated interface.

7.3 Dynamic Alignment

Coronal Plane

Ensure that M-L thrust is minimal by adjusting relative positions of the socket and foot.

Sagittal Plane

Check for smooth transition from heel strike to toe-off. Ensure that when standing the heel and toe are evenly loaded and that both are touching the floor.

⚠ Check that the serrations are engaged before tightening the foot bolt. Incorrectly engaged serrations can lead to noise and result in sudden loosening of the foot. Confirm correct engagement of serrations by slackening foot bolt 1/8 turn. Foot should be loose but not rotate, then tighten to correct torque.

Important: ensure the device is orientated correctly with the clamping screw positioned posterior.

* Allow for users own footwear

8 Fitting Advice

The instructions in this section are for practitioner use only.

934331 for snubber replacement

| Symptoms | Cause | Solution |
|---|-------------------|--|
| Sinking at heel strike Foot slap | Heel too soft | Confirm A-P alignment and heel height accommodation Fit a harder ankle ball |
| Rapid transition from the heel strike through stance phase Difficulty controlling heel action (e.g. foot jars into mid-stance) Foot feels too rigid | Heel too hard | Confirm A-P alignment and heel height accommodation Fit a softer ankle ball |
| Drop off in mid-late stance Feels like walking downhill | Forefoot too soft | Confirm A-P alignment and heel height accommodation Fit a harder snubber |
| Difficulty walking over the forefoot Feels like walking uphill | Forefoot too hard | Confirm A-P alignment and heel height accommodation Fit a softer snubber |

Note... Bilateral amputees may require stiffer snubbers than unilateral amputees. Trans-femoral amputees may require a softer plantarflexion action than trans-tibial amputees.

Ankle Ball and Snubber Selection



Snubbers are sold separately.

| User Weight | 60 kg (131 lb) | 80 kg (176 lb) | 100 kg (220 lb) | 125 kg (275 lb) |
|---------------------|-------------------|-------------------|--------------------|--------------------|
| Ankle Ball | | | | |
| Ankle Part Number | 409314 | 409314 | 409315 | 409316 |
| Snubber | | | | |
| Snubber Part Number | 409007 | 409008 | 409008 | 409009 |

W = White

O = Orange


B = Blue

9 Technical Data

| | |
|--|--|
| Operating and Storage Temperature Range: | -10°C to 50°C (14°F to 122°F) |
| Component Weight: | 199 g (7 oz) |
| Activity Level: | 1–2 |
| Maximum User Weight: | Levels 1–3: 125 kg (275 lb) Level 4: 100 g (220 lb) |
| Attachment Type: | Proximal—Blatchford Male Pyramid Adapter Distal—Serrated end with 3/8" threaded hole for Blatchford Multiflex foot bolt |
| Range of Adjustment: | ±7° tilt from vertical |
| Build Height: | 58 mm |

Build Height

10 Ordering Information

| Part | Part Number |
|----------------------------|-------------|
| Multiflex Ankle W | 409353 |
| Multiflex Ankle O | 409354 |
| Multiflex Ankle B | 409355 |
| Snubber 60° W | 409007 |
| Snubber 70° O | 409008 |
| Snubber 80° B | 409009 |
| Ankle ball joint W | 409314 |
| Ankle ball joint O | 409315 |
| Ankle ball joint B | 409316 |
| Collet | 400702 |
| Circlip | 918000 |
| Rubber and nylon lubricant | 928107 |

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

The product should be recycled where possible in accordance with local waste recycling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Multiflex and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.