Epirus

Instructions for Use

EP22L1S-EP30R8S
EP22L1SD-EP30R8SD

EN	Instructions for Use	2
SL	Navodila za uporabo	17
BG	Инструкции за употреба	32
HR	Upute za uporabu	47
SK	Návod na použitie	62
HU	Használati útmutató	77
EL	Οδηγίες χρήσης	92
LV	Lietošanas pamācība	107
LT	Naudojimo instrukcija	122
ET	Kasutusjuhend	137

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	5
3 Construction	6
4 Function.....	7
5 Maintenance	7
6 Limitations on Use	8
7 Bench Alignment.....	9
7.1 Static Alignment.....	9
7.2 Dynamic Alignment	9
8 Fitting Advice.....	10
9 Assembly Instructions	11
9.1 Fitting/Replacing Buffer and Stop Pad.....	11
9.2 Foot Shell Removal.....	12
9.3 Device and Foot Shell Assembly.....	12
10 Technical Data	14
11 Ordering Information	15

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to Epirus.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for single user.

A moderate-energy-return foot with resilient multi-axial ankle movement. Independent heel and toe spring provides some axial deflection. The split toe provides good ground compliance.

Activity Level

This device is recommended for users that have the potential to achieve Activity Level 3 (weight limits apply, see *Technical Data* section).

Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances. There may also be a number of users in the Activity Levels 2 and 4* who would benefit from the smoother transition from heel to toe offered by this device, but this decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

*Maximum user weight 100 kg and always use one higher spring rate category than shown in the Spring set selection table.

Contraindications

This device might not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Clinical Benefits

- Improved energy return from springs
- Improved mediolateral ground compliance

The remaining instructions in this section are for practitioner use only.

Spring Set Selection

This is a guide for initial spring & plantar flexion buffer selection, refer to Fitting Advice section for final choice.

Activity Level 3

44-52 (100-115)	53-59 (116-130)	60-68 (131-150)	69-77 (151-170)	78-88 (171-195)	89-100 (196-220)	101-116 (221-255)	117-125 (256-275)	kg (lbs)	User Weight
1	2	3	4	5	6	7	8		Foot Spring set
None		Medium			Firm				Suggested Plantar Flexion Buffer

Notes:

- If in doubt choosing between two categories, choose the higher rate spring set.
- Foot Spring set recommendations shown are for transtibial users.
- For transfemoral users or lower activity users we suggest selecting a spring set one category lower, refer to Section 8 *Fitting Advice* to ensure satisfactory function and range of movement.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately. Ensure any use of the device complies with the conditions given in *Limitations on Use*.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.

Avoid exposure to extreme heat and/or cold.

The user must not adjust or tamper with the setup of the device.

The user should contact their practitioner if their condition changes.


Be aware of finger trap hazard at all times.

3 Construction

Principal Parts

- Carrier Assembly Aluminum/stainless steel
- Pyramid Stainless steel
- Pyramid Screw Titanium
- Dome Aluminum
- Plantar Flexion Buffer Natural rubber
- Spring Attachment Screws Titanium/stainless steel
- Heel & Toe Springs e-carbon
- Glide Sock UHM PE (Ultra-High-Molecular-Weight Polyethylene)
- Foot Shell PU (Polyurethane)

Component Identification

4 Function

This device comprises an e-carbon toe and independent heel spring combined with resilient multi-axial ankle movement, ankle function can be customised to the individuals requirements via the use of interchangeable buffers.

Heel and toe springs are attached to the carrier using titanium and stainless steel screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or restricted/excessive rotation, significant wear or excessive discoloration from long term exposure to UV.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Remove the foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness (see *Construction* section), clean and reassemble as necessary.
- Visually check the heel and toe springs for signs of delamination or wear and replace if necessary. Some surface damage may occur after a period of use, this does not affect the function or strength of the foot.
- Remove the dome, check the condition of the plantar flexion buffer and stop pad, replace if necessary (see *Fitting/Replacing Buffer and Stop Pad* section).

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the foot is recommended and any changes in performance or signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing this device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.


Suitable for outdoor use

Exclusively for use between -15 °C and 50 °C
(5 °F to 122 °F)..

7 Bench Alignment

The instructions in this section are for practitioner use only.

7.1 Static Alignment

Set-up Length

With flexion, adduction and abduction properly accommodated, set the limb length 5 mm longer than the sound side to allow for compression and deflection of the foot springs during gait. This should be re-assessed once the dynamic trial commences and the length adjusted accordingly.

Build Line

This should fall between the centre line of the pyramid and 10 mm anterior (with heel height properly accommodated). The socket should be positioned accordingly.

7.2 Dynamic Alignment

Coronal Plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

Sagittal Plane

Check for smooth transition from heel strike to toe-off. Ensure also that when standing the heel and toe are **evenly loaded** and that both are **touching** the floor.

8 Fitting Advice

The instructions in this section are for practitioner use only.

Spring sets are supplied as matched pairs i.e. the heel and toe spring are designed to work together. Progression throughout the stance phase should be smooth; the combined heel and ankle function is key to this process.

If a smooth progression is not achieved, the following should be considered:

	Symptoms	Remedy
Ankle too soft	<ul style="list-style-type: none">• Rapid plantar flexion and possible foot slap• Lack of energy return from heel spring	Fit plantar flexion buffer. A choice of medium and firm buffer is available (the medium version is castellated), see figs 1 & 2 below. For trial purposes the buffer does not need to be bonded. Final selection should be bonded in position using Loctite 424.
Heel spring too soft	<ul style="list-style-type: none">• Sinking at heel strike• Difficulty climbing over the toe• (toe feels too hard)	<ol style="list-style-type: none">1. Move socket anteriorly in relation to the foot (excess movement may result in drop-off)2. Fit a stiffer spring set
Heel spring too hard	<ul style="list-style-type: none">• Rapid transition from heel strike through stance phase• Difficulty in controlling heel action, foot jars into mid-stance• Foot feels too rigid	<ol style="list-style-type: none">1. Move socket posteriorly in relation to foot2. Fit softer spring set

Please contact your supplier if it is not possible to achieve a smooth gait after following the advice above.

9 Assembly Instructions

The instructions in this section are for practitioner use only.

 Use appropriate health and safety equipment at all times including extraction facilities.

 Be aware of finger trap hazard at all times.

9.1 Fitting/Replacing Buffer and Stop Pad

9.1.1 Fitting/Replacing Buffer

1. Unscrew pyramid screw and remove pyramid screw, pyramid and dome.
2. Remove completely any existing buffer from 'C'-shaped recess on underside of dome.
3. Degrease mating surfaces of the 'C'-shaped recess and new buffer.
4. Apply Loctite 424 adhesive to the mating surfaces and bond new buffer into the recess.
5. When the bond has set, refit dome, pyramid and pyramid screw.
6. Tighten pyramid screw to 20Nm.

Fig. 1. Dome Removed Ready for Fitting Buffer



Fig. 2. Medium Buffer Fitted

9.1.2 Fitting/Replacing Stop Pad




1. Unscrew pyramid screw and remove pyramid screw, pyramid and dome.
2. Remove completely any worn stop pad from the square recess on underside of dome.
3. Lightly abrade and degrease mating surfaces of the the new stop pad and square recess.
4. Apply Loctite 424 adhesive to the mating surfaces and bond new stop pad into the recess.
5. When the bond has set, refit dome, pyramid and pyramid screw, and apply a smear of PTFE-loaded grease between stop pad and carrier assembly.
6. Tighten pyramid screw to 20Nm.

9.2 Foot Shell Removal

 **Be aware of finger trap hazard at all times.**

9.3 Device and Foot Shell Assembly

5

If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.

6

Slide carrier/heel spring assembly into the foot shell.

7

Toe spring location in foot shell.

8

Use a suitable lever to encourage the heel spring into the location in the foot shell.

9

Ensure heel spring is engaged into slot.

Heel spring location slot

10


Ensure glide sock does not get trapped when assembling to female pyramid part.

10 Technical Data

Operating and Storage Temperature Range:	-15°C to 50°C (5°F to 122°F)
Component Weight (Size 26N):	610 g (1 lb 6 oz)
Activity Level:	2–3
Maximum User Weight:	Levels 1–3: 125 kg (275 lb) Level 4: 100 kg (220 lb)
Proximal Alignment Attachment:	Male Pyramid (Blatchford)
Range of Adjustment	+/- 7° angular
Build Height: (See diagram below)	85 mm
Heel Height	10 mm

Fitting Length

Size	A
22–26	65 mm
27–28	70 mm
29–30	75 mm

11 Ordering Information

Order Example

EP	25	L	N	3	S
	Size	Side (L/R)	Width* (N/W)	Spring Set Category	Sandal Toe

Available from size 22 to size 30:

EP22L1S to EP30R8S

EP22L1SD to EP30R8SD

(add 'D' for a dark tone foot shell)

*Sizes 25-28 only. For all other sizes, omit the Width field.

e.g. EP25LN3S, EP22R4S, EP27RW4SD

Spring Kits				
Rate	Foot Sizes			
	Small (S)	Medium (M)	Large (L)	Extra Large (XL)
	22-24	25-26	27-28	29-30
Set 1	539801S	539810S	539819S	539828S
Set 2	539802S	539811S	539820S	539829S
Set 3	539803S	539812S	539821S	539830S
Set 4	539804S	539813S	539822S	539831S
Set 5	539805S	539814S	539823S	539832S
Set 6	539806S	539815S	539824S	539833S
Set 7	539807S	539816S	539825S	539834S
Set 8	539808S	539817S	539826S	539835S

Foot Shell (for dark add 'D')		
Size/Side	Narrow	Wide
22L	539038S	-
22R	539039S	-
23L	539040S	-
23R	539041S	-
24L	539042S	-
24R	539043S	-
25L	539044SN	539044SW
25R	539045SN	539045SW
26L	539046SN	539046SW
26R	539047SN	539047SW
27L	539048SN	539048SW
27R	539049SN	539049SW
28L	539050SN	539050SW
28R	539051SN	539051SW
29L	-	539052S
29R	-	539053S
30L	-	539054S
30R	-	539055S

Item	Part. No.
Glide Sock (Sizes 22-26)	531011
Glide Sock (Sizes 27-30)	532811
Plantar Flexion Buffer: Medium	533710
Plantar Flexion Buffer: Firm	533711
Stop Pad	533708

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

The device is warranted for 36 months - foot shell 12 months - glide sock 3 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Epirus and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

 Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH.

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Blatchford Europe GmbH
Endolite India Ltd.