Ø 30 mm Transtibial Long Stump Adapter Kit

Instructions for Use

339240

EN	Instructions for Use	2
SL	Navodila za uporabo	12
BG	Инструкции за употреба	22
HR	Upute za uporabu	32
SK	Návod na použitie	42
HU	Használati útmutató	52
EL	Οδηγίες χρήσης	62
LV	Lietošanas pamācība	72
LT	Naudojimo instrukcija	82
ET	Kasutusjuhend	92

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	6
5 Maintenance	6
6 Limitations on Use	7
7 Fitting Advice.....	8
8 Technical Data	9
9 Ordering Information	10

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to the Ø 30mm Transtibial Long Stump Adapter Kit.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

The device is designed for use exclusively as part of a lower limb prosthesis.

The device is approved for use with Blatchford modular components and thermoplastic or laminated sockets.

Constructed in accordance with Blatchford socket manufacturing recommendations.

Intended for single user only.

Features

When tightened securely the device offers a secure way to connect two parts of the limb. The limb system in conjunction with other components allows for angular, rotational and length adjustments.

Activity Level

This device is suitable for activity levels 1 to 4; weight limits apply, see Technical Data, Section 8. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

There are no known contraindications if used in accordance with these instructions.

Clinical Benefit

- Allows for a low build height connection between prosthetic componentry.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Only combine this device with Modular components and Thermoplastic or Laminated Sockets.

Any changes in the performance or function of the device e.g. unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Avoid exposure to extreme heat and/or cold.

To reduce the risk of injury due to failure or loosening of the screw connections ensure the bolt threads are cleaned thoroughly before each installation.

Be aware of finger trap hazard at all times.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician that has attended an approved training course.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their practitioner if their condition changes.

Always apply the specified torque value to bolt. Never replace the bolt with an alternative bolt. Always use the bolt specified.


Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

3 Construction

Principal Parts

- Coupling Assembly Aluminum
- Attachment Bolt Stainless steel
- Friction Disc Synthetic rubber
- Socket Plate Stainless Steel
- Long Stump Adaptor 30mm Tube Aluminum
- Eccentric Washer Stainless Steel

Component Identification

4 Function

This device is designed to securely connect two parts of the limb build.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider.

Changes in performance may include:

- Any unusual noise
- Excessive play or loss alignment

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

Maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

Carry out the following routine maintenance at least annually:

- Ensure all screws are secure. If not remove and clean screws, reapply Loctite and tighten all screws to the correct torque settings.
- Check for visual defects that may affect proper function.

6 Limitations on Use

This device must only be used with Blatchford modular components and thermoplastic or laminated sockets with Blatchford spherical end form.

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

The product is waterproof to a maximum depth of 1 metre.

This device is suitable for submersion in fresh water only. Light surface corrosion affects neither the function nor the security of this device. However, if it is heavily corroded, stop using and contact your practitioner.

Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to the device.

Thoroughly rinse with fresh water after use in salt or chlorinated water.

Exclusively for use between -15 °C and 50 °C
(5 °F to 122 °F).



Suitable for submersion

7 Fitting Advice

The instructions in this section are for practitioner use only.

1. Cut Adapter square to required length using Cutting Jig 941237. Remove any burrs from the cut end before assembly.
2. Loosely assemble to socket and mark surplus bolt length protruding from socket plate to ensure full engagement when tightened.
3. Cut bolt to length, apply Loctite 222 (926001) to threads and reassemble. Torque tighten to 35 Nm.


Symptom	Remedy
A recurring noise occurs at the interface between the socket and this device.	Ensure that the socket has been constructed in accordance with Blatchford recommendations. Ensure that the socket isn't failing at the connection interface with this device. Confirm the security, Loctite and torque of the device. (Refer to 934366 - Instruction manual for laminated sockets). Ensure the tube is cut square to fit.
The adapter moves out of position.	User must not use the device until adjusted, repaired or replaced.

8 Technical Data

Operating and Storage Temperature Range:	-15°C to 50°C (5°F to 122°F)
Component Weight (<i>size 26</i>):	185 g (7 oz)
Activity Level:	1–4
Maximum User Weight:	100kg (220lb) 1–4 125kg (275lb) 1–3
Attachment Type:	Proximal—Single Bolt Alignment Distal—Ø 30mm Tube
Range of Adjustment:	360° rotation ±6° angular ±4 mm shift
Build Height:	See diagram below

Build Height

Storage and Handling

When storing for prolonged periods, ensure the product is free from moisture and stored at room temperature.

9 Ordering Information

Part Number

339240

Spare Parts	
BK Upper Alignment Coupling	409090
BK Long Stump Adapter 30mm diameter	334003
Washer 30mm Long Stump Adapter	334004
Foot Bolt 3/8" UNF	390102
Socket Plate SYMES	390103
Peg Spanner	940060

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

The above may vary by market; consult your local representative for details.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is constructed from recyclable materials. Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.