Compact SAKL

Instructions for Use

P239143

| | | |
|----|-----------------------|----|
| EN | Instructions for Use | 2 |
| DA | Brugsanvisning | 15 |
| NO | Bruksanvisning | 28 |
| FI | Käyttöohjeet | 41 |
| SV | Bruksanvisning | 54 |
| RU | Инструкция протезиста | 67 |
| ZH | 使用说明 | 80 |
| AR | إرشادات الاستخدام | 93 |

Contents

| | |
|--|----|
| Contents | 2 |
| 1 Description and Intended Purpose | 3 |
| 2 Safety Information | 4 |
| 3 Construction | 5 |
| 4 Function..... | 5 |
| 5 Maintenance | 6 |
| 6 Limitations on Use | 7 |
| 7 Bench Alignment..... | 8 |
| 8 Fitting Instructions..... | 9 |
| 8.1 Unlocking Procedure | 10 |
| 9 Assembly Instructions | 11 |
| 10 Technical Data | 12 |
| 11 Ordering Information | 13 |

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to Compact Semi-Automatic Knee Lock (Compact SAKL).

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

The device is a lightweight knee, ideally suited to individuals with low activity level who need the security of a locked knee during walking.

Features

- Slim, low profile knee
- 140° knee flexion
- Locks with an audible 'click' on full extension
- Thigh release

Activity Level

This device is recommended for users that have the potential to achieve Activity Level 1 or 2. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Contraindications

Not for use with a turntable above the device.

Clinical Benefits

- Lightweight device provides stable ambulation.
- Automatically locks in extension for security when walking.
- For sitting, thigh release cable provides an easy and accessible method of releasing knee into flexion.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their practitioner if their condition changes.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk.

Be aware of finger trap hazard at all times


Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

3 Construction

Principal Parts

- Pyramid Stainless steel
- Chassis Aluminum
- Pivot Pin Stainless steel
- Body Assembly Aluminum, hardened steel, phosphor bronze
- Screw and Washer Stainless steel
- Patella Retainer Polyurethane
- Patella Nylon

Component Identification

4 Function

The device locks in full extension with an audible click allowing secure walking with a locked knee.

The device remains locked in full extension until the release lever is depressed, allowing the knee to bend for sitting purposes. The release lever, attached to the prosthetic socket, operates the release mechanism of the device through a connecting cable.

Note... The knee lock should not be released whilst walking.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or restricted/excessive rotation, significant wear or excessive discoloration from long term exposure to UV.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check all screws for tightness, clean and reassemble as necessary
- If ML rock is evident, the blue plastic bearings 923700 may be renewed, see *Assembly Instructions*
- If AP rock is evident, the buffers 234310 may be renewed, see *Assembly Instructions*
- If the cable is stretched, remove the operating lever knob and loosen the cable clamp set screw.
With the lever up, take up slack in cable without operating the pawl and retighten the set screw after applying Loctite 222. Cut excess cable if necessary, leaving enough for future adjustment, fold back into lever and refit knob.
See also 938120, *Thigh Release* instructions for use.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the foot is recommended and signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

Any changes in performance of this device must be reported to the practitioner

Changes in performance may include:

- Reduced support during stance
- Any unusual noise or movement
- Stiffness or failure to release from lock in thigh release

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits. Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).


Suitable for outdoor use

7 Bench Alignment

The instructions in this section are for practitioner use only.

For socket manufacture use ESK height position on Dummy Knee Block on alignment mandrel.

Alignment adjustments can affect cable tension and consequently knee release.





Always check the knee locks and functions as intended following alignment adjustments.

8 Fitting Instructions

The instructions in this section are for practitioner use only.

 Use appropriate health and safety equipment at all times including extraction facilities.




 Be aware of finger trap hazard at all times.

8.1 Unlocking Procedure

The instructions in this section are for practitioner and user.

When releasing the lock ensure a flexion load is not being applied and the foot is on the ground.

9 Assembly Instructions

The instructions in this section are for practitioner use only.

1 Bearing Replacement

234313

Loctite 243

STP

6

10Nm

234312

923700

923700

2 Patella Replacement

CLICK!

P234314

234315

PULL

3 Buffer Replacement

Do NOT Remove patella from pivot.

Patella Retainer

Buffer 234310

Loctite 424

Note... Patella retainer band may require replacement if damaged during buffer replacement process.

4 Cable Replacement

Release pivots to replace cable if necessary
(See 1 Bearing Replacement)

239243


1

2

10 Technical Data

| | |
|--|-----------------------------------|
| Operating and Storage Temperature Range: | -15 °C to 50 °C 5 °F to 122 °F |
| Component Weight: | 325 g (12 oz) |
| Recommended Activity Level: | 1/2 |
| Maximum User Weight: | 125 kg (275 lb) |
| Maximum Knee Flexion Angle: | 140° |
| Proximal Angular Adjustment: | ±7° |
| Distal Angular Rotation: | 360° |
| Proximal Alignment attachment: | Male Pyramid (Blatchford) |
| Distal attachment: | 30 mm Diameter |
| Build Height: | 40 mm |
| Proximal pyramid to distal tube end | |
| [See diagram below] | |

Fitting Length

11 Ordering Information

Order Example

| Part | Part Number |
|--------------|-------------|
| Compact SAKL | P239143 |

Replacement Parts

| Part | Part Number |
|---------|------------------|
| 923700 | Plastic Bearing |
| P234314 | Patella |
| 234315 | Patella Retainer |
| 234310 | Buffer |

Accessories

| Part | Part Number |
|--|-------------|
| Thigh Release Kit (Release Lever Only) | P239443 |
| Thigh Release Cable | 239243 |

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.