

Smart IP

Blatchford was the first company in the world to have a microprocessor knee (MPK) available on the prosthetic market; the Intelligent Prosthesis (IP) in 1993. Sensors were incorporated and took measurements to determine walking speed and the behaviour of the pneumatic section of the piston was adapted, providing the appropriate degree of swing phase extension. A further iteration of this technology included the IP+, which had programming advancements and was simpler for the prosthetist to calibrate. The algorithm changed to no longer simply measure time between full extensions of the limb, but also used knee flexion information too.

Smart IP is an intelligent prosthetic knee with mechanical control in stance (via an ESK - most frequently supplied with the stance flexion version) and microprocessor-controlled pneumatic swing control. It builds on the previous IP knee iterations, with the same mechanical structure, however it introduced a new learning algorithm, allowing the prosthesis to self-adjust settings as the user walks.

Improvements in Clinical Outcomes using prosthetic knees with microprocessor-controlled swing phase

Improvement in **SAFETY**

- Less cognitive demand during walking, leading to reduced postural sway¹

Improvement in **MOBILITY**

- Increased walking speed²⁻⁵
- Easier to walk at different speeds^{4,6}
- More natural gait⁴
- Easier to walk on slopes^{4,6}

Improvement in **ENERGY EXPENDITURE**

- Reduced energy expenditure compared to (non-MPK) mechanical knees³⁻⁸
- Equivalent energy expenditure to other MPKs (swing *and* stance controlled)⁹
- Reduced self-perceived effort^{4,6}
- Energy expenditure closer to that of able-bodied control subjects¹⁰
- Able to walk further before becoming tired⁴

Improvement in **SYMMETRY**

- Better step length symmetry^{2,6}

Improvement in **USER SATISFACTION**

- Preference over other prosthetic knees^{4,6}

References

1. Heller BW, Datta D, Howitt J. A pilot study comparing the cognitive demand of walking for transfemoral amputees using the Intelligent Prosthesis with that using

- conventionally damped knees. *Clin Rehabil* 2000; 14: 518–522.
2. Chin T, Maeda Y, Sawamura S, et al. Successful prosthetic fitting of elderly trans-femoral amputees with Intelligent Prosthesis (IP): a clinical pilot study. *Prosthet Orthot Int* 2007; 31: 271–276.
 3. Datta D, Heller B, Howitt J. A comparative evaluation of oxygen consumption and gait pattern in amputees using Intelligent Prostheses and conventionally damped knee swing-phase control. *Clin Rehabil* 2005; 19: 398–403.
 4. Datta D, Howitt J. Conventional versus microchip controlled pneumatic swing phase control for trans-femoral amputees: user's verdict. *Prosthet Orthot Int* 1998; 22: 129–135.
 5. Buckley JG, Spence WD, Solomonidis SE. Energy cost of walking: comparison of "intelligent prosthesis" with conventional mechanism. *Arch Phys Med Rehabil* 1997; 78: 330–333.
 6. Kirker S, Keymer S, Talbot J, et al. An assessment of the intelligent knee prosthesis. *Clin Rehabil* 1996; 10: 267–273.
 7. Chin T, Sawamura S, Shiba R, et al. Energy expenditure during walking in amputees after disarticulation of the hip: a microprocessor-controlled swing-phase control knee versus a mechanical-controlled stance-phase control knee. *J Bone Joint Surg Br* 2005; 87: 117–119.
 8. Taylor MB, Clark E, Offord EA, et al. A comparison of energy expenditure by a high level trans-femoral amputee using the Intelligent Prosthesis and conventionally damped prosthetic limbs. *Prosthet Orthot Int* 1996; 20: 116–121.
 9. Chin T, Machida K, Sawamura S, et al. Comparison of different microprocessor controlled knee joints on the energy consumption during walking in trans-femoral amputees: intelligent knee prosthesis (IP) versus C-leg. *Prosthet Orthot Int* 2006; 30: 73–80.
 10. Chin T, Sawamura S, Shiba R, et al. Effect of an Intelligent Prosthesis (IP) on the walking ability of young transfemoral amputees: comparison of IP users with able-bodied people. *Am J Phys Med Rehabil* 2003; 82: 447–451.