



Paediatric Footwear

Ordering Information

1. Ordering and Enquiries:

Ordering Paediatric footwear from Blatchford is simple and easy. You can contact our Customer Services directly:

Telephone: +44 (0) 114 2637900

Email: orthoticsales@blatchford.co.uk

Facsimile: +44 (0) 114 2637901

Post: Blatchford Limited, 11 Atlas Way, Atlas North, Sheffield, S4 7QQ.

2. Delivery Promise:

Our goal is to provide fast and efficient delivery. We aim to despatch all orders for Paediatric footwear from receipt in 14 working days. This does not apply where specific delivery schedules have been previously agreed.

3. Catalogues/Product Information:

If you require additional catalogues or product information please contact our customer service advisors who will be happy to help with your enquiry.

4. Shortages/Damaged Goods:

If any items are received damaged or there are shortages in the delivery this must be notified in writing to customer services within 5 days of receipt of goods.

5. Returns:

In the unlikely event of any footwear being supplied faulty or having been incorrectly supplied, please contact customer services so that the correct or replacement goods may be despatched immediately and a collection arranged for the items to be returned. We are happy to accept Paediatric Footwear that is returned for credit if they are not suitable for the patient.





















For items to be accepted for credit the following conditions apply:

- The footwear must be unworn and returned in less than 30 days from the date of the delivery
- The product must not have been adapted in any way
- Split sized made to order footwear is non-returnable for credit
- Split sized footwear direct from stock is returnable for credit
- Are returned in the original packaging
- Are in a re-saleable condition

Blatchford reserves the right to reject items returned for credit if the above conditions are not met. Blatchford is unable to accept return of items that are 'made to order' for credit, unless the goods are faulty or it can be demonstrated they have been damaged in transit.


Products and prices are subject to variation at any time and are subject to availability. E & OE.

Contents

	Page
Mendivil Paediatric Footwear	3
<hr/>	
Stocked Footwear	4
 Slickers Boot II	4
 Catapult Boot II	4
 Sprinter Boot	4
 Cub Starter Boot	4
<hr/>	
Stock to Order Footwear	5
Specification Guide/Adaptations and Insoles	5-6
 Shikra	7
 Owl	7
 Raptor	8
 Eagle	8
 Kite	9
 Osprey	9
 Condor	10
 Harrier	10
 Merlin	11
 Peregrine	11
 Slickers Boot II	12
 Catapult Boot II	12
 Sprinter Boot	13
 Cub Starter Boot	13
<hr/>	
Orthoses Accommodating Footwear	14
 Falcon	14
 Hawk	14
<hr/>	
Mendivil A-Line Boot	15-16
<hr/>	
Measurement Charts	17-19
Ordering Information	20

Colour Options

Not all colour options are applicable to every style, please check the specification provided alongside each boot in the catalogue pages.



Leather - Box
Colour - White



Leather - Nubuck
Colour - Ademuz Arce



Leather - Nubuck
Colour - Lago



Leather - Nubuck
Colour - Orchid



Leather - Natural
Colour - Black



Leather - Nubuck
Colour - Pink



Leather - Nubuck
Colour - Tartufo



Leather - Bolonia
Colour - Red



Leather - Nubuck
Colour - Chine



Leather - Nubuck
Colour - Black



Leather - Bolonia
Colour - Pink




Leather - Nubuck
Colour - Light Brown



Leather - Bolonia
Colour - Lilac



Leather - Natural
Colour - Royal



Leather - Patent
Colour - Black




Leather - Nubuck
Colour - Fuchsia




Leather - Val
Colour - Soan



Leather - Calgary
Colour - Lilac



Leather - Nubuck
Colour - Navy



Leather - Calgary
Colour - Black




Leather - Bolonia
Colour - Fuchsia



Leather - Val
Colour - Saddale



Leather - Nubuck
Colour - Brown



Leather - Antic
Colour - Jean

The four Stocked Footwear styles on this page are available for 48 hour delivery in the option colours and fastening shown.

Cub Starter Boot



Sizes 18-26
Last Medium
Fastening Lay On Hook and Loop
Height UK
Stiffeners Standard
Colour Nubuck Navy

PU Sole  Black

Slickers Boot II



Sizes 27-45
Last Medium and Wide
Fastening Lace
Height UK
Stiffeners Standard
Colour Natural Leather Black

Rubber Montana  Black

Catapult Boot II



Sizes 27-45
Last Medium and Wide
Fastening Lace
Height UK
Stiffeners Standard
Colour Nubuck Black

Rubber Montana  Black

Sprinter Boot




Sizes 27-45
Last Medium and Wide
Fastening Lace
Height UK
Stiffeners Standard
Colour Nubuck Black

Sports  Black

Stock to Order Specification Guide

Please clearly state the following options when ordering your chosen style. Not all options are applicable to every style, please check the specification provided alongside each boot in the catalogue pages.

- 1 Choose from either a **lace** or **hook and loop** fastening (hook and loop incurs an additional charge).
- 2 Specify the correct height:



- 3 **5 standard last** and **3 orthoses accommodating last** options are available. See page 18 for measurement charts. (Split sizes are available but will incur an additional charge.)

- 4 **3 stiffener** options are available:



- 5 **5 sole options** are available:



- 6 **Various colour and material options** can be specified for each style, as shown here.

- 7 All standard width MTO footwear is supplied with the following internal inlay.



Adaptations and Insoles

A17-PAEDS	High boot uppers additional 25mm to footwear
B1	Outside compensating raise up to 25mm with a rand
B2	Each subsequent 25mm or part thereof to B1
B3	Outside compensating raise up to 25mm without a rand
B4	Each subsequent 25mm or part thereof to B3
B5	Raising the heel in a single item of footwear up to 25mm
B7	Addition of a heel seat or valgus/varus wedge fitted to shoe (each)
B8-PAEDS	Stirrup, heel or caliper retaining strap (each)
B9	Any additional padding in a single item of footwear
B15-PAEDS	Addition of orthosis stiffener
B16	Addition of a metatarsal bar fitted to the sole (each)
B17	Inclusion of a rocker sole (each)
B18	Inclusion of wedge, fitted into/onto the sole (each)
B19	Inclusion of a floated sole in a single item of footwear
B20	Inclusion of a sole strip or plate in composite material (each)
B21	Inclusion of a sole strip or plate in non-composite material (each)
B22	Addition of a rubber tip to the sole (each)
B23	Addition of a tip, with flange, to protect the welt (each)
B24	Addition of an external protective toe cap (each)
B25	Full/half-elongated, oblique, floated heel or heel wedge or torque heel in a single item of footwear
B30	Socket to one/both sides of the heel, without front/backstop to single item of footwear
B31	Socket to one/both sides of the heel, with front/backstop to single item of footwear
B32	Inclusion of a toe raising socket to single item of footwear
B34-PAEDS	Addition of a lined and padded double-decker T or Y strap (each)
B35	Rebuilding the heel to allow socket fixation in a single item of footwear
B40a	3/4 or full length non-moulded inlay, non-shock absorb material e.g. Leather/Neoprene/Plastazote® (each)
B40b	3/4 or full length non-moulded inlay, shock absorb material e.g. Poron®/PPT® (each)
B40c	3/4 or full length non-moulded inlay in silicone type material (each)
B41	Metatarsal/valgus or heel pad to inlay (each)
B42	Metatarsal/valgus or heel pad fitted direct to shoe (each)
B43	Inclusion of a heel elevator in cork, HDP or EVA covered (each)
B50	Filler to compensate for toe amputation(s) (each)
B51	Filler to compensate for amputation of forepart of foot (each)
NSB030	Flat carbon fibre insole (each)
WLN-PAEDS	Addition of lambswool lining in a single item of footwear (each) (The dark blue lines on the diagram below illustrate the area lined.)



Shikra (Product Code 711)




Last Options 5 standard widths and 3 orthoses widths
Sizes 20-44
Fastening Lace or Reverse Pull Hook and Loop (with medial zip)
Height UK
Stiffeners Standard, Stability/High Boot and Orthosis Accommodation

PU Sole




Black White Pink

Rubber Kick



Black Light Beige



Leather - Bolonia Colour - Pink



Leather - Natural Colour - Royal



Leather - Bolonia Colour - Lilac



Leather - Bolonia Colour - Red



Leather - Nubuck Colour - Orchid



Leather - Nubuck Colour - Tartufo

Owl (Product Code 712)



Last Options 5 standard widths and 3 orthoses widths
Sizes 16-45
Fastening Lace or Reverse Pull Hook and Loop (with medial zip)
Height UK
Stiffeners Standard, Stability/High Boot and Orthosis Accommodation

PU Sole



Black Grey Pink



Leather - Nubuck Colour - Fuchsia



Leather - Nubuck Colour - Black



Leather - Nubuck Colour - Chine



Leather - Nubuck Colour - Tartufo



Leather - Natural Colour - Black



Leather - Natural Colour - Royal

Raptor (Product Code 713)



Last Options 5 standard widths and 3 orthoses widths
Sizes 18-45
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Stability/High Boot and Orthosis Accommodation



Leather - Calgary
Colour - Lilac

Leather - Calgary
Colour - Black



Leather - Bolonia
Colour - Lilac

Leather - Natural
Colour - Black



Leather - Nubuck
Colour - Tartufo

Leather - Natural
Colour - Royal

Eagle (Product Code 715)



Last Options 5 standard widths and 3 orthoses widths
Sizes 20-44
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Stability/High Boot and Orthosis Accommodation



Leather - Calgary
Colour - Lilac

Leather - Calgary
Colour - Black



Leather - Nubuck
Colour - Pink Upper/
Ademuz Arce Cuff

Leather - Bolonia
Colour - Fuchsia



Leather - Nubuck
Colour - Light Brown

Leather - Nubuck
Colour - Fuchsia

Kite (Product Code 716)



Last Options 5 standard widths and 3 orthoses widths
Sizes 20-44
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Stability/High Boot and Orthosis Accommodation



Leather - Bolonia
 Colour - Red Upper/
 White Stitching



Leather - Nubuck
 Colour - Fuchsia Upper/
 White Stitching



Leather - Natural
 Colour - Black Upper/
 Cream Stitching



Leather - Nubuck
 Colour - Black Upper/
 Cream Stitching



Leather - Natural
 Colour - Royal Upper/
 White Stitching



Leather - Val
 Colour - Saddale/
 Brown Stitching


Osprey (Product Code 717)



Last Options 5 standard widths and 3 orthoses widths
Sizes 20-44
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Stability/High Boot and Orthosis Accommodation




Leather - Nubuck
 Colour - Fuchsia Upper/
 Ademuz Arce Cuff



Leather - Patent
 Colour - Black/
 Ademuz Arce Cuff



Leather - Calgary
 Colour - Black



Leather - Calgary
 Colour - Lilac



Leather - Bolonia
 Colour - Lilac Upper/
 Fuchsia Nubuck Cuff



Leather - Natural
 Colour - Royal

Condor (Product Code 718)




- Last Options** 5 standard widths and 3 orthoses widths
- Sizes** 18-45
- Fastening** Lace or Reverse Pull Hook and Loop (with medial zip)
- Height** Continental
- Stiffeners** Standard, Stability/High Boot and Orthosis Accommodation

Rubber Montana



Leather - Calgary
Colour - Black



Leather - Calgary
Colour - Lilac



Leather - Patent
Colour - Black



Leather - Bolonia
Colour - Lilac



Leather - Natural
Colour - Black



Leather - Natural
Colour - Royal

Harrier (Product Code 719)



- Last Options** 5 standard widths and 3 orthoses widths
- Sizes** 18-45
- Fastening** Lace or Reverse Pull Hook and Loop
- Height** UK or Continental
- Stiffeners** Standard, Stability/High Boot and Orthosis Accommodation

PU Sole




Rubber Montana



Leather - Nubuck
Colour - Chine Upper/
Nubuck Lago Cuff



Leather - Patent
Colour - Black



Leather - Calgary
Colour - Black



Leather - Calgary
Colour - Lilac



Leather - Natural
Colour - Black




Leather - Natural
Colour - Royal

Merlin (Product Code 724)



Last Options 5 standard widths and 3 orthoses widths
Sizes 20-45
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Stability/High Boot and Orthosis Accommodation



Peregrine (Product Code 725)



Last Options 5 standard widths and 3 orthoses widths
Sizes 20-45
Fastening Lace or Reverse Pull Hook and Loop
Height UK or Continental
Stiffeners Standard, Stability/High Boot and Orthosis Accommodation



Slickers Boot II (Product Code 78000)



Last Options	5 standard widths
Sizes	18-45
Fastening	Lace or Reverse Pull Hook and Loop
Height	UK
Stiffeners	Standard, Stability/High Boot and Orthosis Accommodation



Leather - Natural
Colour - Black



Leather - Val
Colour - Saddale



Leather - Nubuck
Colour - Black



Leather - Val
Colour - Soan



Leather - Natural
Colour - Royal



Leather - Nubuck
Colour - Pink

Catapult Boot II (Product Code 78020)



Last Options	5 standard widths
Sizes	18-45
Fastening	Lace or Reverse Pull Hook and Loop
Height	UK
Stiffeners	Standard, Stability/High Boot and Orthosis Accommodation



Leather - Natural
Colour - Black



Leather - Val
Colour - Saddale



Leather - Nubuck
Colour - Black



Leather - Natural
Colour - Royal



Leather - Nubuck
Colour - Tartufo



Leather - Nubuck
Colour - Pink

Sprinter Boot (Product Code 78500)



Last Options 5 standard widths
Sizes 20-45
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Stability/High Boot and Orthosis Accommodation

Sports



Cub Starter Boot (Product Code 70706)




Last Options 5 standard widths
Sizes 16-26
Fastening Lace or Reverse Pull Hook and Loop
Height UK
Stiffeners Standard, Stability/High Boot and Orthosis Accommodation

PU Sole



Leather - Nubuck
Colour - Black



Leather - Box
Colour - White



Leather - Val
Colour - Saddale



Leather - Natural
Colour - Black



Leather - Natural
Colour - Royal



Leather - Nubuck
Colour - Pink



Leather - Nubuck
Colour - Black



Leather - Nubuck
Colour - Navy



Leather - Nubuck
Colour - Fuchsia



Leather - Nubuck
Colour - Pink



Leather - Val
Colour - Soan



Leather - Nubuck
Colour - Orchid

Orthoses Accommodating Footwear

The styles below feature an additional zip fastening for ease of access and donning.

Falcon (Product Code 714)



- Last Options** 3 orthoses widths
- Sizes** 20-44
- Fastening** Lace with Zip
- Height** UK
- Stiffeners** Standard, Stability/High Boot and Orthosis Accommodation

Rubber Montana  Black (with yellow stitch effect)

Rubber Kick  Black (with yellow stitch effect)



Leather - Calgary
Colour - Black



Leather - Calgary
Colour - Lilac



Leather - Natural
Colour - Black



Leather - Nubuck
Colour - Black



Leather - Nubuck
Colour - Brown




Leather - Natural
Colour - Royal

Hawk (Product Code 726)



- Last Options** 3 orthoses widths
- Sizes** 20-45
- Fastening** Lace or Reverse Pull Hook and Loop with Zip
- Height** UK
- Stiffeners** Standard, Stability/High Boot and Orthosis Accommodation

Sports    Black Grey White



Leather - Box
Colour - White Upper/
Grey Cross Detail



Leather - Natural
Colour - Black Upper/
Grey Cross Detail




Leather - Natural
Colour - Royal Upper/
Grey Cross Detail



Leather - Bolonia
Colour - Red



Leather - Bolonia
Colour - Fuchsia



Leather - Bolonia
Colour - Pink



Mendivil A-Line Boot

- A-Line boot designed for the effective treatment of Congenital Talipes Equino Varus (CTEV)
- Manufactured using softer leather giving greater comfort, tolerance and compliance
- Reversible padded tongue
- Increased calcaneal cupping and padding
- Posterior cut out facilitating accurate positioning and monitoring of the heel
- Elongated fastening strap with clearly stitched markers
- Standard and narrow widths available
- Works in conjunction with a click bar using the Ponseti method



Model No. 900
Standard or
905 Narrow

Product Specifications

Please specify required fitting and size.

Mendivil A-Line Boot Standard Model No. 900 and 907-S	Colour White	Sizes 12-33
Mendivil A-Line Boot Narrow Model No. 905	Colour White	Sizes 12-33
Mendivil A-Line Boot Narrow Model No. 907-N	Colour White	Sizes 12-33

Model No. 907
Narrow or
Standard



Mendivil A-Line Boot Standard - Model No. 900 Mendivil A-Line Boot Tall - Model No. 907-S

Size	12	13	14	15	16	17	18	19	20	21	22	23	24
Joint Width	46.00	47.00	48.00	50.00	51.00	52.00	53.00	55.00	56.00	57.00	59.00	60.00	61.00
Length	73.00	79.00	85.00	92.00	99.00	105.00	112.00	119.00	125.00	133.00	139.00	146.00	154.00
Joint Circumference	128.00	131.50	135.00	138.50	142.00	145.50	149.00	152.50	156.00	159.50	163.00	166.50	170.00
Size	25	26	27	28	29	30	31	32	33				
Joint Width	62.00	63.00	64.00	66.00	67.00	68.00	70.00	71.00	72.00				
Length	159.00	167.00	173.00	180.00	187.00	193.00	201.00	208.00	215.00				
Joint Circumference	173.50	177.00	180.50	184.00	187.50	191.00	194.50	198.00	201.50				

Mendivil A-Line Boot Narrow - Model No. 905 Mendivil A-Line Boot Tall - Model No. 907-N

Size	12	13	14	15	16	17	18	19	20	21	22	23	24
Joint Width	40.00	42.00	43.00	45.00	46.00	48.00	49.00	50.00	52.00	54.00	55.00	56.00	58.00
Length	73.00	79.00	85.00	92.00	99.00	105.00	112.00	119.00	125.00	133.00	139.00	146.00	154.00
Joint Circumference	116.00	119.50	123.00	126.50	130.00	133.50	137.00	140.50	144.00	147.50	151.00	154.50	158.00
Size	25	26	27	28	29	30	31	32	33				
Joint Width	59.00	61.00	62.00	63.00	65.00	67.00	68.00	69.00	70.00				
Length	159.00	167.00	173.00	180.00	187.00	193.00	201.00	208.00	215.00				
Joint Circumference	161.50	165.00	168.50	172.00	175.50	179.00	182.50	186.00	189.50				



Mendivil Paediatric Footwear

Mendivil Paediatric Footwear features a range of styles in sizes 16 to 45 with two availability options: Stocked Footwear or Stock to Order Footwear. Stocked Footwear is available for a 48 hour delivery and limited to the four styles shown on page 3. These four styles are offered in all sizes and in standard or wide width fittings, but only in the colour and fastening options shown.

Made to Order Footwear styles are shown from page 4 onwards. This range offers a wide selection of colour options in either lace or hook and loop fastening in two back heights; UK and Continental. Standard and orthoses accommodating width fittings plus three stiffening options are also available. Not all options are applicable to every style, please check the specification provided alongside each boot in the catalogue on the relevant page.

All Mendivil Paediatric Footwear can be fully adapted with a range of floats, wedges and raises to meet your clinical requirements.



Measurement Charts

Fitting Principles

The measurement charts lay out a step by step procedure which can help the clinical assessment and measurement process achieve a right first time fit.

It is important to note that when the foot weight bears, the foot posture changes and can have an effect on the shoe size required. Please consider this to ensure the desired toe room is added to the overall length.

1. Foot length
2. Joint width
3. Joint circumference
4. Instep circumference

Orthosis Accommodating Last

Boot		Orthoses Narrow Last			Orthoses Wide Last			Orthoses Extra Volume		
Size	Length	Joint Width	Joint Circ.	Long Heel	Joint Width	Joint Circ.	Long Heel	Joint Width	Joint Circ.	Long Heel
20	133	57	170	225	60	178	235	60	184	241
21	139	59	174	233	62	182	243	62	188	249
22	146	60	178	241	63	186	251	63	192	257
23	153	61	182	249	64	190	259	64	196	265
24	159	62	186	257	65	194	267	65	200	273
25	166	64	190	265	67	198	275	67	204	281
26	173	66	194	273	69	202	283	69	208	289
27	179	68	199	281	71	207	291	71	213	297
28	186	70	203	289	73	211	299	73	217	305
29	193	71	207	297	74	215	307	74	221	313
30	199	72	212	305	75	220	315	75	226	321
31	206	73	217	315	76	225	325	76	231	331
32	213	75	221	323	78	229	333	78	235	339
33	219	77	225	331	80	233	341	80	239	347
34	226	78	229	339	81	237	349	81	243	355
35	232	80	233	347	83	241	357	83	247	363
36	239	82	237	355	85	245	365	85	251	371
37	246	83	240	365	86	248	375	86	254	381
38	252	85	244	375	88	252	385	88	258	391
39	259	86	248	385	89	256	395	89	262	401
40	266	88	252	395	91	260	405	91	266	411
41	272	90	255	405	93	263	415	93	269	421
42	279	91	258	415	94	266	425	94	272	431
43	286	93	261	425	96	269	435	96	275	441
44	292	95	264	435	98	272	445	98	278	451
45	299	96	268	444	99	275	454	99	281	460

All measurements are in mm for Length, Joint Width, Instep Circumference and Long Heel Circumference.



Standard Last

Boot		Narrow		Medium		Wide		X Wide		XX Wide	
Size	Length	Joint Width	Instep Circ.	Joint Width	Instep Circ.	Joint Width	Instep Circ.	Joint Width	Instep Circ.	Joint Width	Instep Circ.
16	106	46	122	49	130	52	138	55	146	58	154
17	113	47	126	50	134	53	142	56	150	59	158
18	119	48	130	51	138	54	146	57	154	60	162
19	126	50	134	53	142	56	150	59	158	62	166
20	133	51	138	54	146	57	154	60	162	63	170
21	139	53	142	56	150	59	158	62	166	65	174
22	146	54	146	57	154	60	162	63	170	66	178
23	153	55	150	58	158	61	166	64	174	67	182
24	159	56	154	59	162	62	170	65	178	68	186
25	166	58	158	61	166	64	174	67	182	70	190
26	173	60	162	63	170	66	178	69	186	72	194
27	179	62	166	65	174	68	182	71	190	74	198
28	186	64	170	67	178	70	186	73	194	76	202
29	193	65	174	68	182	71	190	74	198	77	206
30	199	66	178	69	186	72	194	75	202	78	210
31	206	67	182	70	190	73	198	76	206	79	214
32	213	69	186	72	194	75	202	78	210	81	218
33	219	71	190	74	198	77	206	80	214	83	222
34	226	72	194	75	202	78	210	81	218	84	226
35	232	74	198	77	206	80	214	83	222	86	230
36	239	76	202	79	210	82	218	85	226	88	234
37	246	77	206	80	214	83	222	86	230	89	238
38	252	79	210	82	218	85	226	88	234	91	242
39	259	80	214	83	222	86	230	89	238	92	246
40	266	82	218	85	226	88	234	91	242	94	250
41	272	84	222	87	230	90	238	93	246	96	254
42	279	85	226	88	234	91	242	94	250	97	258
43	286	87	230	90	238	93	246	96	254	99	262
44	292	89	234	92	242	95	250	98	258	101	266
45	299	90	238	93	246	96	254	99	262	102	270

All measurements are in mm for Length, Joint Width, Instep Circumference and Long Heel Circumference.

Distributed by

Blatchford:

+44 (0) 114 2637900 | orthoticsales@blatchford.co.uk
Blatchford Limited, 30 Atlas Way, Sheffield, S4 7QQ, United Kingdom.

@blatchfordgrp | blatchford.co.uk

224183998 Iss1 12/19. Information correct at time of print.

