

AIRCAST®

XP Diabetic Walker™

The XP Diabetic Walker System™ is specifically designed to meet the needs of the higher risk diabetic patient, providing protection and immobilization. It makes an ideal alternative to total contact casting, eliminating the need for costly and time consuming removal and application during skin and wound examinations. The XP Diabetic Walker can be easily removed and reapplied time after time and closed securely with the Aircast Brace-Lok™ to help monitor patient progress. A specialized Hand Bulb with a pressure gauge is included to measure individual aircell compression, ensuring regulated inflation for neuropathic patients.

- Multi-aircell support for a 'total contact' fit
- Rocker sole for reduction of plantar pressures and improved off-loading
- Dual density insole to help eliminate pressure points
- Graduated aircell compression for edema control
- Regulated aircell inflation using the specialized Hand Bulb with pressure gauge
- Semi-rigid shell combined with full aircell coverage maximizes plantar unloading, regulates shear stress, and provides protection and immobilization
- Extension straps are available to increase the strap length by 4" (10 cm)



	Men's Shoe Size	Women's Shoe Size	Height of Brace	Maximum Foot Width at Toes	Maximum Calf Circumference
Small	4-7	5-8	32cm (12.5")	11cm (4.31")	36cm (14")
Medium	7-10	8-11	35cm (13.75")	12cm (4.688")	41cm (16")
Large	10-13	11-15	38cm (14.875")	13cm (5.25")	66cm (26")

Technical support is available from Steve Osborne, please email steveo@blatchford.co.uk

blatchford.co.uk

Chas A Blatchford & Sons Ltd., 11 Atlas Way, Atlas North, Sheffield, S4 7QQ.

Tel: +44 (0) 114 2637900 Fax: +44 (0) 114 2637901 Email: orthoticsales@blatchford.co.uk Web: www.blatchford.co.uk

