

Hallux-Valgus Correctors



Orliman Nighttime Hallux-Valgus Corrector



Product Codes: NSE846 (Right) & NSE847 (Left)

- Encircles the entire foot
- Velcro fastening provides ease of adjustment
- Incorporates a pliable aluminium abductor splint that can be shaped to regulate the position of the big toe
- Comfortable cotton bouclé lining

Orliman Daytime Hallux-Valgus Corrector



Product Codes: NSE924

- Utilises a cap on the big toe that is held securely in place with adjustable elastic banding around the heel
- Velcro fastening on the banding allows for control over traction to achieve accurate positioning
- Supportive elastic banding runs around the instep

blatchford.co.uk