

- Provides optimal cervical support for both children and adults
- Non toxic hypoallergenic plastazote foam
- Range of heights and circumferences available from stock - see chart below
- Velcro fasteners
- Completely radiolucent
- Easily cleaned with soap and water
- Optional snap on Halo and Stabiliser attachments for extra immobilisation are available to order*



Wheaton Philadelphia Collar

Product Code: NSM053

Size	Circumference	Height 2"	Height 2¼"	Height 3¼"	Height 4¼"	Height 5¼"
Infant	6-8"	✓	✗	✗	✗	✗
Paediatric	8-11"	✗	✓	✗	✗	✗
Small	10-13"	✗	✓	✓	✓	✓
Medium	13-16"	✗	✓	✓	✓	✓
Large	16-19"	✗	✓	✓	✓	✓
Extra Large	19"+	✗	✓	✓	✓	✓

*Please request Halo and Stabiliser as required additions when ordering.

To order please contact your Customer Services Representative or alternatively, email orthoticsales@blatchford.co.uk

blatchford.co.uk

Chas A Blatchford & Sons Ltd., 11 Atlas Way, Atlas North, Sheffield, S4 7QQ.

Tel: +44 (0) 114 2637900 Fax: +44 (0) 114 2637901 Email: orthoticsales@blatchford.co.uk Web: www.blatchford.co.uk