


VASYLI+
Dananberg



The very latest in 1st ray function orthotic technology

The Vasyli + Dananberg insole is the innovative result of a collaboration between Vasyli and Howard J. Dananberg, named amongst the most influential podiatrists in America by Podiatry Management Magazine. It features removable proximal and distal plugs that can be detached when the orthotic is fitted to coincide with the head of the 1st MPJ. In cases where maximum effect is required, both plugs can be removed, offering a larger path of least resistance to the joint in an effort to enhance 1st ray function. A selection of additions are also available to further complement the orthotic.



- 1st ray technology with removable proximal and distal plugs
- Suede finish cover impregnated with Antibacterial FootFresh®
- Enhanced metatarsal control
- Heat mouldable
- Features a unique Tri-Compound Construction (TCC™). This blend of PU, EVA and Sorbon ensures comfort, function and control

Product Code	Product Size	UK Sizes		European Sizes	
		Men	Women	Men	Women
71HD04	Extra Small (XS)	3-4½	3½-5	35½-38½	35½-38½
71HD06	Small (S)	5-6½	5½-7	39-40	39-40
71HD08	Medium (M)	7-8½	7½-9	40½-42½	40½-42½
71HD10	Large (L)	9-10½	9½-11	43-45	43-45
71HD12	Extra Large (XL)	11-12½	-	46-48	-