



How to Order

To order any product featured in this leaflet, please contact your Customer Services Representative or alternatively, email orthoticsales@blatchford.co.uk



The Orliman catalogue, featuring a complete range of orthoses is available upon request. Please contact John Stephen on 0114 263 7910 or email johnst@blatchford.co.uk

Knee Orthoses

blatchford.co.uk

Chas A Blatchford & Sons Ltd., 11 Atlas Way, Atlas North, Sheffield, S4 7QQ.
Tel: +44 (0) 114 2637900 Fax: +44 (0) 114 2637901 Email: orthoticsales@blatchford.co.uk Web: www.blatchford.co.uk

Cover montage incorporates images © Chrisarvey & © Sebastian Kaultzki | Agency: Dreamstime.com 960029 0908



Three Panel Knee Immobiliser

Product Codes: NSE820 (50cm)/NSE822 (60cm)

- Immobilises the knee, retaining the required position
- Three panel system adjusts to accommodate the individual wearer
- Fastening straps on the interior panel provide a secure fit
- Beige velour fabric exterior with cotton bouclé lining
- Universal size available in 50cm and 60cm length



Knee Orthosis with Flexion/Extension Joint & Airlock Block

Product Code: NSE963

- Monocentric jointed knee orthosis that provides stability and control over a range of motion
- Comprises two femoral and two tibial adjustable foam cuffs with velcro straps
- Block in extension at 0°, or for the control of flexion/extension from 0° to 80° in 10° intervals
- Available in black
- Universal size available in 43cm and 59cm length



Adjustable Knee Orthosis with Lock

Product Code: 94260



Easily adjustable and lockable without the need for tools!

Exceptional value!

- Fully adjustable knee orthosis with monocentric articulation and locking system
- Provides good functionality, stability and control over a range of motion
- Easily adjustable joint - no tools required
- Easily extendable from 44cm to 69cm, using telescopic bars and friction clips
- Comprises two femoral and two tibial cuffs with velcro straps
- 0° to 120° flexion control in 10° intervals
- -10° to 30° extension control in 10° intervals
- Padded mediolateral plates for improved comfort
- Available in black
- Universal size

