

Sizing Charts

Please note: there is no exact translation between UK and European sizing. The charts below represent the closest size equivalent.



78010/1/2 AFC Boot II

Narrow

Size (Euro/UK)	27/9½	28/10	29/11	30/12	31/12½	32/13½	33/1	34/2	35/2½	36/3½
Joint Circ.	176.00	180.20	184.40	188.60	192.80	197.00	201.20	205.40	209.60	213.80
Joint Width	65.00	66.60	68.20	69.80	71.40	73.00	74.60	76.20	77.80	79.40
Toe Depth	21.00	21.60	22.20	22.80	23.40	24.00	24.60	25.20	25.80	26.40

Size (Euro/UK)	37/4	38/5	39/6	40/6½	41/7½	42/8	43/9	44/10	45/10½
Joint Circ.	218.00	222.20	226.40	230.60	234.80	239.00	243.20	247.40	251.60
Joint Width	81.00	82.60	84.20	85.80	87.40	89.00	90.60	92.20	93.80
Toe Depth	27.00	27.60	28.20	28.80	29.40	30.00	30.60	31.20	31.80

Standard

Size (Euro/UK)	27/9½	28/10	29/11	30/12	31/12½	32/13½	33/1	34/2	35/2½	36/3½
Joint Circ.	184.00	188.20	192.40	196.60	200.80	205.00	209.20	213.40	217.60	221.80
Joint Width	67.00	68.60	70.20	71.80	73.40	75.00	76.60	78.20	79.80	81.40
Toe Depth	21.50	22.10	22.70	23.30	23.90	24.50	25.10	25.70	26.30	26.90

Size (Euro/UK)	37/4	38/5	39/6	40/6½	41/7½	42/8	43/9	44/10	45/10½
Joint Circ.	226.00	230.20	234.40	238.60	242.80	247.00	251.20	255.40	259.60
Joint Width	83.00	84.60	86.20	87.80	89.40	91.00	92.60	94.20	95.80
Toe Depth	27.50	28.10	28.70	29.30	29.90	30.50	31.10	31.70	32.30

Wide

Size (Euro/UK)	27/9½	28/10	29/11	30/12	31/12½	32/13½	33/1	34/2	35/2½	36/3½
Joint Circ.	192.00	196.20	200.40	204.60	208.80	213.00	217.20	221.40	225.60	229.80
Joint Width	69.00	70.60	72.20	73.80	75.40	77.00	78.60	80.20	81.80	83.40
Toe Depth	22.00	22.60	23.20	23.80	24.40	25.00	25.60	26.20	26.80	27.40

Size (Euro/UK)	37/4	38/5	39/6	40/6½	41/7½	42/8	43/9	44/10	45/10½
Joint Circ.	234.00	238.20	242.40	246.60	250.80	255.00	259.20	263.40	267.60
Joint Width	85.00	86.60	88.20	89.80	91.40	93.00	94.60	96.20	97.80
Toe Depth	28.00	28.60	29.20	29.80	30.40	31.00	31.60	32.20	32.80



78017 AFC Over Splint Boot

Medium

Size (Euro/UK)	27/9½	28/10	29/11	30/12	31/12½	32/13½	33/1	34/2	35/2½	36/3½
Joint Circ.	197.00	201.30	205.60	209.90	214.20	218.50	222.80	227.10	231.40	235.70
Joint Width	68.00	69.50	71.00	72.50	74.00	75.50	77.00	78.50	80.00	81.50
Toe Depth	26.00	26.80	27.60	28.40	29.20	30.00	30.80	31.60	32.40	33.20

Size (Euro/UK)	37/4	38/5	39/6	40/6½	41/7½	42/8	43/9	44/10	45/10½
Joint Circ.	240.00	244.30	248.60	252.90	257.20	261.50	265.80	270.10	274.40
Joint Width	83.00	84.50	86.00	87.50	89.00	90.50	92.00	93.50	95.00
Toe Depth	34.00	34.80	35.60	36.40	37.20	38.00	38.80	39.60	40.40

Choice and Delivery

AFC Boots in Navy Blue leather with a grain effect finish, as illustrated in this brochure, can be delivered from stock on a next working day basis. They can also be made to order in a choice of colour and leather options at no additional cost and delivered within 10 to 15 working days from receipt of the customer order. All AFC Boots can be supplied with split sizes and on a 'mix and match' basis e.g. left foot, medium fitting, size 29 combined with right foot, wide fitting, size 30. Please see our current Orthotic Products Price List for full terms and conditions.

For more information on the AFC Boot or any other products or service from Blatchford Blatchford please contact Customer Services. Alternatively, contact John Stephen on +44 (0) 114 263 7910 or email johnst@blatchford.co.uk




Paediatric AFC Boot

...for the world's most precious feet





Blue Nubuck



Red Nubuck



Black Nubuck



Tan Leather



Red Leather



Navy Blue Leather



Pink Leather

Colour Options Please note that due to print limitations the actual leather colours may vary slightly from the samples below.

AFC Boot II



AFC Over Splint Boot



Paediatric AFC Boot

Mendivil Blatchford AFC Boots are supplied with the practical addition of a functional foot orthosis as standard, providing a clinically effective footwear solution that offers exceptional value for money. The AFC Boot has recently been refined and re-styled to offer even better performance and value. The introduction of an Over Splint AFC Boot gives extra depth where required and an extended selection of colour and leather options further enhances the range. A variety of adaptations and repairs are also available at an additional cost.

The AFC Boot II is manufactured using our new range of specifically designed paediatric shoe lasts. The foot length measurements are comparable to other leading footwear brands.

- Three width fittings available
- Low throat opening for ease of donning and doffing
- High quality materials
- Hard wearing sole unit
- Padded collar and tongue
- AFC Orthoses included as standard
- Exceptional value for money

AFC Orthoses

AFC Boot functional foot orthoses are heat mouldable to afford the best fit, superior comfort and greater compliance.

- Three quarter medium density EVA shell (shore 60)
- Three quarter low density EVA shell (shore 40) available upon request
- Medium density EVA full length top cover
- Integrated 3° varus rear foot posting
- Additional postings available (to specification)

More complex functional foot orthoses can be produced utilising Exocam state of the art CAD/CAM manufacturing technology. Please contact us for more details.

