


Insole Top Cover Options



Anti-odour, Thermodynamic and Anti-microbial

Blatchford offers a diverse range of top cover options on all custom made insoles:

- Poron (grey, green & pink)
- PPT (backed & unbacked)
- X2 memory
- Spenco
- Silano microfibre (gold & black)
- Cambrelle
- Leather/Perforated Leather
- Vinyl
- Nora
- Plastazote
- Black sponge rubber
- X-static



X-static® is available as a top cover option on all Blatchford custom made insoles. It is an innovative fibre that has a natural layer of pure silver permanently bonded to it. The silver is irreversibly bound to a polymer so it becomes a physical part of the fibre rather than acting as a surface treatment. Silver provides unparalleled anti-odour, thermodynamic, and anti-microbial properties that will not diminish over time. More than 10 years of research and testing stand behind the superb performance of X-Static®, it is safe, natural and will last for the life of the product.

To order, please contact your Customer Services Representative or alternatively, email orthoticsales@blatchford.co.uk

blatchford.co.uk