

Wheaton Pavlik Harness™

Specially designed to allow control over hip movement without the rigidity and discomfort of casting.



- Velcro fasteners allow for rapid, yet secure adjustment
- Colour coding on anterior straps helps parents to avoid misapplication
- Anti-slip foot pieces are specially designed to prevent the foot from slipping out of the harness
- The harness has no obtrusive buckles and the halter has a soft inner lining so the child can be comfortable in any position
- Carefully made to conform to prescribed orthopaedic standards
- All adjustments including the posterior straps can be made while the child is supine

For single patient use only

blatchford.co.uk