

# Becker Lerman Haro 2000™

*Hip Abduction Rotation Orthosis indicated following dislocation after total hip replacement or hip revision surgery.*



- Anterior pelvic section facilitates ease of donning, comfort is enhanced by eliminating pressure on the back and spine
- Integral adhesion layers on the pelvic and thigh liners are designed to control rotation and prevent migration of the orthosis
- Innovative attachment of joint to pelvic and thigh modules allows rotational adjustment and ensures accurate anatomical alignment
- Allows adjustment and accurate alignment in transverse plane
- Offers control and adjustment of flexion/extension and abduction/adduction

Model Numbers	Size	Waist
145 LS/145 RS	Small	28-32
145 LM/145 RM	Medium	34-38
145 LL/145RXL	Large	40-44
145 LXL/145 RXL	Extra Large	46-50

**blatchford.co.uk**