

Vasyli Sizing Chart

Product	Product Code	Intrinsic Rearfoot Varus Posting	Green Sizes Low Density	Grey Sizes Low/Medium Density	Blue Sizes Medium Density	Red Sizes High Density
Footprint Regular 3/4 Length	NSB300	4°	XS to XL	n/a	XS to XL	Kids to XL
Footprint Full Length	NSB310	4°	XS to XL	n/a	XS to XL	Kids to XL
Footprint Slimfit	NSB320	2°	n/a	n/a	XS to M	n/a
Footprint Easyfit	NSB300	4°	n/a	n/a	Kids to XL	n/a
Custom Regular 3/4 Length	NSB300	6°	Kids to XL	n/a	S Kids to XL	S Kids to XL
Custom Full Length	NSB310	6°	Kids to XL	n/a	S Kids to XL	S Kids to XL
Custom Slimfit	NSB320	2°	n/a	n/a	XS to M	n/a
Custom Easyfit	NSB300	6°	n/a	n/a	Kids to XL	n/a
Pro Regular 3/4 Length	NSB330	4°	n/a	Kids to XL	n/a	n/a
Pro Full Length	NSB340	4°	n/a	Kids to XL	n/a	n/a
Pro Slimfit	NSB350	2°	n/a	XS to M	n/a	n/a

Please note: there is no exact translation between UK and European sizing. The chart below represents the closest size equivalent.

Vasyli Size References

Size	S Kids	Kids	XS	S	M	L	XL
UK	11-12½	1-2½	3-4½	5-6½	7-8½	9-10½	11-12½
European	29-32	33-35	35½-38½	39-40	41½-42½	43-45	46-48

Vasyli Orthaheel Pain Reliever

Product	Product Code	Sizes
Vasyli Orthaheel Pain Reliever	NSB305	Large (fits male sizes 7½-12 and females 8-12½)
Vasyli Orthaheel Pain Reliever	NSB305	Small (fits male sizes 3½-7 and females 4-7½)

EsCo Sizing Chart

Product	Product Code	Sizes
EsCo Paediatric Insoles	NSE931	20-42

When ordering insoles, please specify style, density (if applicable), size and product code.

To order from the Blatchford Orthotic Insole range simply call 0114 263 79000 or email orthoticsales@blatchford.co.uk. Alternatively contact your Blatchford Customer Services Representative.



ASYLI Orthotic Insoles

Vasyli Footprint - the Original AOL

- Heat mouldable and compatible with Vasyli Orthotic component additions
- Regular and Full Length available in 3 densities - green (low), blue (medium) and red (high)
- Slimfit and Easyfit available in medium density only
- 4° rearfoot varus posting on Regular, Full Length and Easyfit
- 2° rearfoot varus posting and flat heel area on Slimfit



Slimfit



Easyfit




Regular



Full Length

Vasyli Custom

- Heat mouldable
- Bonded top covers
- Regular and Full Length available in 3 densities - green (low), blue (medium) and red (high)
- Slimfit and Easyfit available in medium density only
- Compatible with Vasyli Orthotic component additions
- 6° rearfoot varus posting and deep heel seat on Regular, Full Length and Easyfit
- 2° rearfoot varus posting and flat heel area with velcro patch for shoe attachment on Slimfit




Regular Red



Full Length Red




Slimfit



Regular Green



Full Length Green



Easyfit



Regular Blue



Full Length Blue

Vasyli Pro


- Non-mouldable
- Bonded top covers
- Low/medium density
- Compatible with Vasyli Orthotic component additions
- 4° rearfoot varus posting and deep heel seat on Regular and Full Length
- 2° rearfoot varus posting and flat heel area on Slimfit



Regular



Full Length



Slimfit

Vasyli Orthaheel Pain Reliever

- Reduces traction and strain on ligaments under the foot giving natural, lasting relief from many types of heel pain
- 4° rearfoot varus angle and deep heel cup
- Shock absorption with a Sorbon Shock Dot in the centre of the heel for added comfort



Vasyli Component Additions

- A full range of self adhesive orthotic additions including:
 - 1 2° and 4° rearfoot varus posts
 - 2 4° and 6° forefoot varus posts
 - 3 4° and 6° forefoot valgus posts
 - 4 Metatarsal raise additions
 - 5 4mm, 6mm and 8mm heel lifts in small, medium or large.

