

Multi-axial ankle function for secure foot placement during gait.

Multiflex Ankle

Ankle Part Number and Description

Snubber Part Number and Description

Max. User Weight/ Activity Level Description

379342 Multiflex ankle W 30mm
 379343 Multiflex ankle O 30mm
 379344 Multiflex ankle B 30mm

409007 - snubber W 60°
 409008 - snubber O 70°
 409009 - snubber B 80°

125kg / 1-3











409353 Pyramid Multiflex W
 409354 Pyramid Multiflex O
 409355 Pyramid Multiflex B

409007 - snubber W 60°
 409008 - snubber O 70°
 409009 - snubber B 80°

125kg / 1-3



Build height	Material	Weight	Ankle and snubber selection
60mm 	Aluminium	201g	<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">  60kg </div> <div style="text-align: center;">  80kg </div> <div style="text-align: center;">  100kg </div> <div style="text-align: center;">  125kg </div> </div> <div style="display: flex; justify-content: space-between; margin-top: 5px;"> <div style="text-align: center;">  Ankle </div> <div style="text-align: center;">  Snubber </div> </div>
58mm 	Aluminium & Stainless Steel	199g	<p>Snubber for standard ankle:</p> 