



Sandal Toe

Product Range

Blatchford 

Foot Range Now With Sandal Toe

Natural look

Flexible footwear choice

Additional sizes*

Available in 2 skin tones



*Additional width sizes are now available for sizes 25, 26 and 27.

Blatchford's new Sandal Toe Footshell offers a more natural look and provides flexibility of footwear choice with its enhanced sandal toe design. It is available in both a Light and Dark shade to help better match skin tones and offers additional width options helping you to find the perfect fit.

Linx

Multi-award winning, integrated limb system, designed to deliver an experience that mimics the incredible and complex structure of the human leg.



Product Code	Size	Side	Width*	Spring Set	Sandal Toe
LINX	25	L	N	3	S

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Echelon^{VAC}

Hydraulic ankle technology paired with an elevated vacuum system for a secure and comfortable socket connection.



Product Code	Size	Side	Width*	Spring Set	Sandal Toe
EVAC	25	L	N	3	S

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Elan

Microprocessor controlled hydraulic ankle for terrain response and speed control.



Product Code	Size	Side	Width*	Spring Set	Sandal Toe
ELAN	25	L	N	3	S

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Echelon

Biomimetic hydraulic ankle control and split toe for stability on uneven or sloping terrain.



Product Code	Size	Side	Width*	Spring Set	Sandal Toe
EC	25	L	N	3	S

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Echelon^{VT}

Biomimetic hydraulic ankle with shock absorber, rotation and split toe for terrain conformance.

Product Code	Size	Side	Width*	Spring Set	Axial Spring	Sandal Toe
ECVT	25	L	N	5	3	S

Example: foot size 25, left, narrow, spring rating 5, axial spring 3, sandal toe.



Javelin

Light weight blade, ideal for moderate impact sports and activities

Product Code	Size	Side	Width*	Spring Set	Sandal Toe	+ Top Housing
JAV	25	L	N	3	S	Adapter

Example: foot size 25, left, narrow, spring rating 3, sandal toe, adapter.



Esprit

Low profile and lightweight foot ideal for longer residual limbs.

Product Code	Size	Side	Width*	Spring Set	Sandal Toe
ESP	25	L	N	3	S

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



Epirus

Multi-axial foot with ankle joint for ground compliance and energy efficient response.

Product Code	Size	Side	Width*	Spring Set	Sandal Toe
EP	25	L	N	3	S

Example: foot size 25, left, narrow, spring rating 3, sandal toe.



*Narrow (N) and Wide (W) available for sizes 25-27 only.
For dark tone add suffix D.

Footshell Width Guide

Use the below table to help choose the correct size when ordering, for sizes 25, 26 and 27 you can now choose between Narrow (N) and Wide (W) fittings.



Size	22	23	24	25	26	27	28	29	30
Footshell Width	76mm	78mm	80mm	82mm (N) 88mm (W)	83mm (N) 91mm (W)	91mm (N) 99mm (W)	100mm	102mm	105mm
Old Footshell	77mm	77mm	80mm	82mm	82mm	87mm	90mm	95mm	101mm

NB: Widths measured at widest point.

Spring Kits and Footshell Information

If you wish to order spring kits or footshells separately, please refer to the order codes below to ensure you receive the correct item.

Spring Kits (for Linx, Elan, Echelon, EchelonVT, EchelonVAC, Esprit, Epirus):

Foot Size	Foot Size	Foot Size	Foot Size
22-24	25-26	27-28	29-30
Set 1 - 539801S	Set 1 - 539810S	Set 1 - 539819S	Set 1 - 539828S
Set 2 - 539802S	Set 2 - 539811S	Set 2 - 539820S	Set 2 - 539829S
Set 3 - 539803S	Set 3 - 539812S	Set 3 - 539821S	Set 3 - 539830S
Set 4 - 539804S	Set 4 - 539813S	Set 4 - 539822S	Set 4 - 539831S
Set 5 - 539805S	Set 5 - 539814S	Set 5 - 539823S	Set 5 - 539832S
Set 6 - 539806S	Set 6 - 539815S	Set 6 - 539824S	Set 6 - 539833S
Set 7 - 539807S	Set 7 - 539816S	Set 7 - 539825S	Set 7 - 539834S
Set 8 - 539808S	Set 8 - 539817S	Set 8 - 539826S	Set 8 - 539835S

(for Javelin Spring Kits):

Foot Size	Foot Size	Foot Size	Foot Size
22-24	25-26	27-28	29-30
Set 1 - 539301S	Set 1 - 539510S	Set 1 - 539519S	Special
Set 2 - 539302S	Set 2 - 539511S	Set 2 - 539520S	Order Only
Set 3 - 539303S	Set 3 - 539512S	Set 3 - 539521S	Set 3 - 539530S
Set 4 - 539304S	Set 4 - 539513S	Set 4 - 539522S	Set 4 - 539531S
Set 5 - 539305S	Set 5 - 539514S	Set 5 - 539523S	Set 5 - 539532S
Set 6 - 539306S	Set 6 - 539515S	Set 6 - 539524S	Set 6 - 539533S
Set 7 - 539307S	Set 7 - 539516S	Set 7 - 539525S	Set 7 - 539534S
Set 8 - 539308S	Set 8 - 539517S	Set 8 - 539526S	Set 8 - 539535S
Set 9 - -	Set 9 - 539518S	Set 9 - 539527S	Set 9 - 539536S

Footshells*:

Foot size Left/Right	Narrow	Wide
22L	539038S	-
22R	539039S	-
23L	539040S	-
23R	539041S	-
24L	539042S	-
24R	539043S	-
25L	539044SN	539044SW
25R	539045SN	539045SW
26L	539046SN	539046SW
26R	539047SN	539047SW
27L	539048SN	539048SW
27R	539049SN	539049SW
28L	-	539050S
28R	-	539051S
29L	-	539052S
29R	-	539053S
30L	-	539054S
30R	-	539055S

*For dark tone add suffix D

+44 (0) 1256 31660 | customer.service@blatchford.co.uk
Blatchford Products Limited, Unit D Antura, Kingsland Business Park,
Basingstoke, RG24 8PZ, United Kingdom.

@blatchfordgrp | blatchford.co.uk

203266220 Iss1 09/19. Information correct at time of print.

Blatchford