


ProTec II Diabetic Insole

The ProTec II Diabetic Insole is a pressure relieving insole designed to decrease the rate of ulceration, enhance comfort and reduce shear force on the foot. It provides the same benefits of the standard ProTec Insole but features the addition of silver bonded fibres that provide improved inhibition of bacterial and fungal growth.



Biomechanical Features

- High memory Poron 96 mid layer
- Compression-set resistance
- Split density EVA base supports stability and balance

Interface Features

- X-STATIC® top cover - antimicrobial fabric clinically proven to reduce HAI-causing contaminants. 99.9% pure metallic silver bonded fibers create an ionic shield, inhibiting the growth of bacteria and fungi. (N.B. not suitable for deep heel cups, for use on shallow trims only)
- Anti-fungal properties
- Thermo regulator

Increased Patient Comfort

- Excellent impact absorption
- Cushioning microcellular structure
- Open cell - breathable

Full insole additions are also available for selection on RMS.

Please choose 'ProTec II Diabetic Insole' in the product section of RMS to specify this insole.