

# ProTib PTTD Insole

The ProTib PTTD Insole is ideal for patients with Tibialis Posterior Dysfunction. Its custom thermoplastic shell gives excellent biomechanical control and the design incorporates a series of shape features to enhance fit and comfort.

## Increased Patient Comfort

- Streamlined trim to provide comfort on heel strike
- Increased comfort rear foot post

## Foot Friendly Shape


- 2mm feather edge for optimum fit
- Outer shell conforms to shoe

## Interface Features

- Durable Spenco Top Cover
- Soft non-abrasive material

## Biomechanical Features

- Medial heel skive and opposing lateral clip
- Reduced profile medial posting
- U shaped 2mm shell profile enhances medial arch strength



Full insole options are available for selection on RMS including skive, heel post and top cover included in the price. Please choose the 'ProTib PTTD Insole' in the product section of RMS to specify this insole.