

Seating
Accessories



Index


Product	Page
Moulded Sandals	3
Shoulder Harness	4
Dynamic Butterfly Harness	5
Single Pull Pelvic Belt	6
Dual Pull Pelvic Belt	6
Pelvic Stabiliser - Single Pull, Metal Buckle	7
Pelvic Stabiliser - Dual Pull, Plastic Buckle	7

BCS Harnesses and belts are not designed to be used as safety devices and must not be used as vehicle restraints for use in transportation. A separate passenger restraint should always be used in line with the recommendation of the wheelchair manufacturer.

Moulded Sandals

Moulded foot supports can assist in achieving the desired seating position or to offer a greater level of control for the feet.

The supports may be located in any position on the wheelchair footrest. Each support is fitted with soft 25mm nylon webbing toe and instep straps as standard, which are easily adjusted and controlled through a quick release buckle.



Size	A	B	C	Product Code
2	80	90	210	NSS903
3	85	100	245	NSS904
4	100	120	280	NSS905

Shoulder Harness

Designed for firmer control to help maintain shoulder and chest alignment.

Available in four sizes, top pull or front pull, with centre width adjustment and lower swivel buckles to keep the pull across the harness consistent and improve comfort.



	Pad Length A	Pad Width B	Pad Width C	Top Pull Product Code	Front Pull Product Code
Small	260	50	40	SHTP - S	SHFP - S
Medium	315	60	50	SHTP - M	SHFP - M
Large	375	70	55	SHTP - L	SHFP - L
X Large	440	80	60	SHTP - XL	SHFP - XL

Dynamic Butterfly Harness

Proven, close fitting support to help maintain postural control and stability.

Manufactured out of a soft neoprene material, this allows sensory feedback from the snug fit, making it feel secure but allowing for some functional movement.

Available in four sizes, along with three female sizes, top pull and front pull, with lower swivel buckles to keep the pull across the harness consistent and improve comfort.




Female	Pad Length A	Pad Width B	Pad Width C	Top Pull Product Code	Front Pull Product Code
Medium	335	37	140	FBTP - M	FBFP - M
Large	395	45	160	FBTP - L	FBFP - L
X Large	440	50	180	FBTP - XL	FBFP - XL

Male	Pad Length A	Pad Width B	Pad Width C	Top Pull Product Code	Front Pull Product Code
Small	250	45	120	BHTP - S	BHFP - S
Medium	325	50	140	BHTP - M	BHFP - M
Large	385	60	160	BHTP - L	BHFP - L
X Large	440	65	180	BHTP - XL	BHFP - XL

Single Pull Pelvic Belt – Metal Buckle

A two part pelvic harness with metal buckle and single centre pull adjustment for flexibility and adjustment.

Manufactured from a soft 9mm neoprene and soft touch webbing for comfort.




	Length A	Length B	Product Code
Small	300	25	PBM - S
Medium	350	38	PBM - M
Large	440	38	PBM - L
X Large	550	38	PBM - XL

Dual Pull Pelvic Belt – Plastic Buckle

A two part pelvic harness with plastic buckle and dual centre pull adjustment for flexibility and a more positive hold on the pelvis.

Manufactured from a soft 9mm neoprene and soft touch webbing for comfort.




	Length A	Length B	Product Code
Small	300	25	PBP - S
Medium	350	38	PBP - M
Large	440	38	PBP - L
X Large	550	38	PBP - XL

Pelvic Stabiliser - Single Pull, Metal Buckle

A two part pelvic stabiliser with metal buckle, single centre pull adjustment for flexibility and two extra straps for extra security. Suitable for clients with moderate to severe positioning needs. The shaping of the pads helps improve the location of the straps for increased comfort.

Manufactured from a soft 9mm neoprene and soft touch webbing for comfort.




	Length A	Length B	Product Code
X Small	260	25	PSM - XS
Small	340	25	PSM - S
Medium	420	38	PSM - M
Large	540	38	PSM - L

Pelvic Stabiliser - Dual Pull, Plastic Buckle

A two part pelvic stabiliser with plastic buckle, dual centre pull adjustment for flexibility and two extra straps for extra security. Suitable for clients with moderate to severe positioning needs. The shaping of the pads helps improve the location of the straps for increased comfort.

Manufactured from a soft 9mm neoprene and soft touch webbing for comfort.



	Length A	Length B	Product Code
X Small	260	25	PSP - XS
Small	340	25	PSP - S
Medium	420	38	PSP - M
Large	540	38	PSP - L

Information

Blatchford posture belt and harnesses are designed to provide postural support and can be used as part of a specialised seating system or standard wheelchair seating.

Safety Information:

Postural harnesses and belts are not designed to be used as safety devices and must not be used as the sole restraint in transportation.

A suitable, separate 3 point passenger seat belt (fastened into a vehicle) must be used to secure the seated occupant during travel in line with the recommendation of the wheelchair manufacturer.

Always ensure that all pelvic belts and harnesses are a snug fit, hold the client securely and are comfortable.

Upper body harnesses should always be fastened around the body. Always allow enough space between the neck and harness to ensure it does not pose a choking risk.

Where an upper body harness is being used, always use in conjunction with a pelvic belt which has been correctly adjusted to minimise the risk of sliding down in the seat, creating a choking risk.

Cleaning:

Blatchford postural belts and harnesses can be cleaned by wiping with a damp cloth and a mild detergent (please do not use bleach or products containing bleach). If using a washing machine we recommend that the temperature is set to 40°C and they are not tumble dried.

Follow Blatchford: | blatchford.co.uk
@BlatchfordGrp | 0114 263 7900
/BlatchfordGrp | seatingsales@blatchford.co.uk

